FlexLink™ TiBase

Instruction For Use - Инструкция за употреба - 使用说明 - Upute za upotrebu - Návod k použití - Brugsanvisning - Gebruiksaanwijzing - Kasutusjuhend - Käyttöohjeet - Mode d'emploi - Gebrauchsleitung - Οδηγίες χρήσης - Használati útmutató - Istruzioni per l'uso - 取扱説明書 - 사용법 - Lietošanas norādījumi - Naudojimo instrukcija - Bruksanvisning - Instrukcja użytkowania - Instruções de uso - Instrucțiuni de utilizare - Инструкция по применению - Návod na použitie - navodila za uporabo - Instrucciones de uso - Bruksanvisning - วิธีการใช้งาน - Kullanma Talimatı - ارشادات الاستخدام

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Disclaimer of liability Elos Medtech Pinol products may only be used according to the instructions and recommendations. The user of Elos Medtech Pinol products should determine their suitability for particular patients and indications. Elos Medtech Pinol disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or damages arising out of or in connection with any errors in professional judgement or practice in the use or placement of the Elos Medtech Pinol products.

Distributed Exclusively by

BIOMET 3i 4555 Riverside Drive
Palm Beach Gardens, FL 33410, USA
1-800-342-5454/ 561-776-6700

Patents: <http://patents.biomet3i.com/patents.pdf>

List of Contents

English - Instruction For Use	3
Български език / Bulgarian - Инструкция за употреба	4
简体中文 / Simplified Chinese - 使用说明.	5
Hrvatski / Croatian - Upute za upotrebu.	6
Čeština / Czech - Návod k použití.	7
Dansk / Danish - Brugsanvisning.	8
Nederlands / Dutch - Gebruiksaanwijzing.	9
Eesti keel / Estonian - Kasutusjuhend.	10
Suomi / Finnish - Käyttöohjeet.	11
Français/ French - Mode d'emploi.	12
Deutsch / German - Gebrauchsanleitung.	13
Ελληνικά / Greek - Οδηγίες χρήσης.	14
Magyar / Hungarian - Használati útmutató	15
Italiano / Italian - Istruzioni per l'uso	16
日本語 / Japanese - 取扱説明書.	17
한국어 / Korean - 사용법.	18
Latviski/ Latvian - Lietošanas norādījumi.	19
Lietuvių k. / Lithuanian - Naudojimo instrukcija.	20
Norsk/ Norwegian - Bruksanvisning	21
Polski/ Polish - Instrukcja użytkowania	22
Português / Portuguese - Instruções de uso (Portugal)	23
Română / Romanian - Instrucțiuni de utilizare	24
Русский / Russian - Инструкция по применению	25
Slovenčina / Slovak - Návod na použitie	26
Slovenščina / Slovenian - navodila za uporabo	27
Español / Spanish - Instrucciones de uso (España)	28
Svenska / Swedish - Bruksanvisning	29
ภาษาไทย / Thai - วิธีการใช้งาน	30
Türkçe / Turkish - Kullanma Talimatı	31
Arشادات الاستخدام العربية / Arab	32

English

Instruction For Use

FlexLink™ TiBase Distributed Exclusively by BIOMET 3i

Indications for Use

The FlexLink™ TiBase is intended for prosthetic restorations. It is used as an interface between a dental implant and an abutment/restoration. The FlexLink™ TiBase will be attached to the implant using a screw and attached to the restoration by cementation. For compatibility FlexLink™ TiBase please refer to information on individual product label.

Device Description

FlexLink™ TiBases are manufactured from biocompatible titanium alloy grade 5 ELI (TiAl6V4 ELI). FlexLink™ TiBases include various surface treatments. For specific product descriptions, please refer to individual product labels.

Contraindications

The FlexLink™ TiBase is not intended for restorations requiring an angle more than 20° to the axis of the implant.

Allergies to the alloy, or contents of the alloy may very rarely occur.

Warnings

FlexLink™ TiBase may only be used for BIOMET 3i implant systems. The connection to the implant system must never be changed or modified. FlexLink™ TiBase and included BIOMET 3i Gold-Tite® or Titanium prosthetic screws are for single use only. Reuse can result in loss of functionality and/or infections. Since the FlexLink™ TiBase and screw are small they must be handled with caution to avoid the patient swallowing or inhaling them.

When mounting the Prosthetic Screw, it is important to use a manual screwdriver before using any kind of torque wrench.

Precautions

The FlexLink™ TiBase must only be used and handled by dental professionals who have approved training in dentistry. The use of abutment screw torque higher than the recommended value may result in damage to the FlexLink™ TiBase, the screw and/or the implant. The use of torque values lower than those recommended may result in loosening of the FlexLink™ TiBase and/or the screw. It is recommended that 3.4mm FlexLink™ TiBase products not be placed in the posterior regions. For any lab work and any try-in of the FlexLink™ TiBase, a lab screw must be used. The included BIOMET 3i Gold-Tite® or Titanium prosthetic screw may only be used to permanently fix the finished restoration.

MRI Statement

The FlexLink™ TiBase has not been evaluated for safety, heating, migration or compatibility in the Magnetic Resonance Imaging (MRI) environment.

Potential Adverse Events

Potential adverse events associated with the use of FlexLink™ TiBase products may include: failure to integrate, loss of integration and infection.

Note for Dental Professionals

Tighten the BIOMET 3i IUNIHG, ILRGHG Gold-Tite® and IUNIHT, ILRGHT, and UNIHT Titanium prosthetic screws to 20Ncm. Tighten the BIOMET 3i UNISG prosthetic screw to 32-35Ncm.

Constructing the Custom Abutment/Restoration

The custom abutment must be constructed using approved design software and FlexLink™ TiBase Scan Posts. Internal geometry of the custom abutment is defined by the library file issued by BIOMET 3i. For a ceramic abutment or restoration, minimum wall thickness is to be met according to the ceramic supplier's recommendation.

Milling a Custom Ceramic Abutment/Restoration

The custom ceramic abutment must be manufactured with milling equipment adequate for the designed geometry of the ceramic abutment.

Sintering the Custom Ceramic Abutment/Restoration

The custom ceramic abutment/restoration should be sintered in accordance with the instruction issued by the ceramic material supplier.

Preparing the Custom Ceramic Abutment/Restoration for Cementing

If the custom ceramic abutment is to be veneered, this must be done in accordance with the instruction issued by the ceramic material supplier. Internal geometry must be blasted in accordance with the instruction issued by ceramic material supplier. Before cementing, the surface should be thoroughly cleaned with alcohol or steam and dried.

Preparing the FlexLink™ TiBase for Cementing

Blasting of FlexLink™ TiBase is not necessary but optional. Before cementing, the surface should be thoroughly cleaned with alcohol or steam and dried. If the FlexLink™ TiBase is to be blasted, it should be done in accordance with the instructions issued by the abutment material supplier. In order to protect the connection to the implant, hold the FlexLink™ TiBase in an implant analog during blasting.Cementing

Follow cementing procedure in accordance with instruction issued by cement supplier. In order to cement the abutment/restoration to the FlexLink™ TiBase, insert the base into the master cast and fix it in place with the waxing screw.

Apply the cement in a circular fashion to the top edge of the contact surface of the FlexLink™ TiBase and to the internal geometry of the abutment/restoration. Rotate the abutment/restoration gently as you push it over the FlexLink™ TiBase until resistance is felt. Then turn it until its final position is reached. The custom abutment/restoration must fit the FlexLink™ TiBase without gaps and voids. Any traces of cement should be removed immediately.

Polishing

Once the cement has set, the excess is removed with rubber polishing tools. The excess cement inside the screw channel must be carefully removed, ideally with the aid of a microscope. When polishing fix the ceramic abutment in an implant analog.

Cleaning and Sterilization

Prior to attachment of the FlexLink™ TiBase/ceramic abutment/restoration and screw to the implant, the FlexLink™ TiBase and Abutment assembly must first be cleaned with water and brush followed by sterilization. The recommended sterilization procedure is full cycle pre-vacuum steam sterilization at a temperature of 135 °C (275°F) for 3 minutes, 4 pulse minimum or full cycle gravity steam sterilization at 135°C (275°F) for 3 minutes minimum. The products are ready for use 5 minutes after autoclaving process.

Caution

US federal law restricts this product for sale by or on order of a dentist or physician. FlexLink™ TiBase abutments may only be used by dentists/ physicians who have had the appropriate education and training.

Ordering Information and Dimensions

See current Product Catalogue for updated and detailed order information of the various FlexLink™ TiBase products.

Storage and Handling

FlexLink TiBase should be stored at room temperature.

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