Elos Abutment Blank

Instruction for use – Gebrauchsanleitung –
Instrucciones de uso – Mode d'emploi –
Istruzioni per l'uso – Gebruiksaanwijzing –
Bruksanvisning

List of content

English – Instruction for use	7
Deutsch – Gebrauchsanleitung	8
Español – Instrucciones de uso	9
Français — Mode d’emploi	10
Italiano – Istruzioni per l’uso.....	11
Nederlands – Gebruiksaanwijzing.....	12
Norsk – Bruksanvisning	14
Svensk - Bruksanvisning.....	15

English

	Catalogue number
	Batch code
	Manufacturer
	Consult instructions for use
	Do not re-use
	Prescription only
	Date of manufacture
	Medical Device

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English – Instruction for use

Indications for use: The Elos Abutment blank is intended for production of patient-specific abutments for attaching to dental implants in order to provide a basis for prosthetic restorations.

Device description: All Elos Abutment Blanks are made from biocompatible Titanium grade 5 (TiAl6V4 ELI).

Contraindication: The Elos Abutment Blank is not intended for restorations requiring an angle more than 30° to the axis of the implant. Allergies to the alloy or contents of the alloy may very rarely occur.

Operating procedure:

Design procedure: The customer-specific abutment must be designed using approved design software. Guidelines for minimum material thicknesses and maximum angles are shown in figure 1. The screw head must not be exposed when designing the customer specific abutment. The handling and use of Elos Abutment Blank requires more than these instructions alone. Only dental surgeons who have undergone approved training in dental implantology should fit the final abutment, only laboratory employees with the relevant training should design the final abutment and only operators with the relevant training should machine the blank for a patient-specific abutment.

Shaping procedure: The customer-specific abutment must be manufactured with milling equipment adequate for the designed geometry. It is important that the product is properly fixed when it is being milled, and also that the Interface of the Abutment Blank is properly protected in order to prevent any damage to it. It is recommended that the Elos Milling Fixture is used to achieve optimal fixing and proper protection.

Clinical procedure: The final abutment is attached to the implant using an Elos Prosthetic Screw and tightened to the recommended torque (refer to individual screw label).
The screw channel must always be sealed after abutment is attached to implant.

Cleaning: Prior to attachment, the final abutment must first be cleaned with water and brush followed by sterilization. The recommended sterilisation procedure is full cycle pre-vacuum steam sterilization at a temperature of 134 °C (273 °F) for 3 minutes minimum or 121°C (250°F) for 15 minutes minimum.

Precautions: In order to ensure the best possible conditions for successful working of an implant, it is strongly recommended that the processing centre that modifies the blank, the laboratory that designs the structure and applies porcelain and the surgeon and dentist who install the components, all work closely together throughout the processing of the implant.

Warnings:

- Only the compatible Elos Prosthetic Screw must be used to secure the Elos Abutment Blank to the implant. Use of other screws than the Elos Prosthetic Screw may result in loss of functionality and/ or infections.
- The use of torque value higher than the recommended may result in damage to the Elos Abutment Blank, the Prosthetic Screw and/or the implant. The use of torque values lower than those recommended may result in loosening of the Elos Abutment Blank, which may result in damage to the Elos Abutment Blank and/or the implant.
- The interface of the Abutment Blank, which being connected to the implant, must never be changed or modified. Modifications of the interface may result in loss of functionality and/ or infections.
- Abutment Blanks are for single use only. Reuse can result in loss of functionality and/ or infections.
- Since the final abutments are small, they must be handled with caution to avoid the patient swallowing or inhaling them. It is recommended that they are handled using forceps


Potential Adverse Events: Potential adverse events associated with the use of Abutment Blank products may include: failure to integrate, loss of integration and infection.

Compatibility: For compatibility of Elos Abutment Blank please refer to information on individual product label.
The Elos Abutment Blank cannot be combined with implants of a different implant type or manufacturer than stated on the label.

Storage and Handling: Elos Abutment Blanks should be stored at room temperature protected from direct sunlight.

Disposal: Must be disposed as biological waste.

MRI Statement: The Elos Abutment Blank has not been evaluated for safety, heating, migration or compatibility in the Magnetic Resonance Imaging (MRI) environment.

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