# Elos Accurate® API® Bridge

Návod k použití – Brugsanvisning –  
Gebrauchsanleitung – Instruction for use –  
Οδηγίες χρήσης – Instrucciones de uso –  
Käyttöohjeet – Mode d'emploi – Használati  
útmutató – Notkunarleiðbeiningar – Istruzioni  
per l'uso – Naudojimo instrukcija –  
Lietošanas norādījumi – Bruksanvisning –  
Gebruiksaanwijzing – Instrukcja użytkowania  
– Instruções de uso – Anvisningar för  
användningen

## List of content

|   |    |
|---|----|
| Čeština - Návod k použití.....              | 3  |
| Dansk - Brugsanvisning.....                 | 5  |
| Deutsch - Gebrauchsanleitung .....          | 7  |
| English - Instruction for use .....         | 9  |
| Ελληνικά - Οδηγίες χρήσης.....              | 11 |
| Español - Instrucciones de uso.....         | 14 |
| Eesti keel - Kasutusjuhised.....            | 16 |
| Suomi - Käyttöohjeet .....                  | 18 |
| Français - Mode d'emploi .....              | 20 |
| Magyar - Használati útmutató .....          | 22 |
| Íslenska - Notkunarleiðbeiningar.....       | 24 |
| Italiano - Istruzioni per l'uso .....       | 26 |
| Lietuvių k. - Naudojimo instrukcija.....    | 28 |
| Latviski - Lietošanas norādījumi.....       | 30 |
| Norsk - Bruksanvisning.....                 | 32 |
| Nederlands - Gebruiksaanwijzing .....       | 34 |
| Polski - Instrukcja użytkowania .....       | 36 |
| Português - Instruções de uso .....         | 39 |
| Svenska - Anvisningar för användningen..... | 41 |

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### Elos Accurate® API® Bridge

#### 1 Použití

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### Elos Accurate® API® Bridge

#### 1 Intended use

The Elos Accurate® API® Bridge is intended for prosthetic restorations. The Hybrid Base™ Bridge is used as an interface between a dental implant or dental abutment and a dental bridge and will be attached to the implant or abutment using the included Prosthetic screw and attached to the dental bridge by cementing.

#### 2 Product Description

The product consists of the Hybrid Base™ Bridge and Prosthetic screw which are manufactured from biocompatible titanium alloy grade 5 ELI (TiAl<sub>6</sub>V<sub>4</sub> ELI). Additionally, it includes a Process screw and a Model analog for preparation of the prosthetic restoration. The product is available for a variety of implant platforms and sizes. For specific product descriptions, please refer to individual product labels.

#### 3 Indications

The Hybrid Base™ Bridge are directly or indirectly connected to the endosseous dental implant and are intended for use as an aid in prosthetic rehabilitations. The dental bridge can be cemented on the Hybrid Base™ Bridge. A temporary restoration can be used prior to the insertion of the final restoration to maintain, stabilize and form the soft tissue during the healing phase. The temporary restoration may not be placed into occlusion. The Hybrid Base™ Bridge may be placed into occlusion when the implant is fully osseointegrated.

#### 4 Contraindications

- The Elos Accurate® API® Bridge is not intended for implants having a divergence angle more than 30° relatively to each other.
- The Elos Accurate® API® Bridge is not intended for single restorations.

#### 5 Warnings and precautions

- The Elos Accurate® API® Bridge components must never be changed or modified.
- The components included in the product are for single-use only.
- Reuse of prosthetic components can result in loss of functionality and/ or infections.
- Tightening and loosening of the included hexalobular Process and Prosthetic screw must be done by use of an Elos Accurate® Prosthetic Screwdriver 18mm, 26mm or 34mm (Ref. No. PS-AH18-1, PS-AH26-1 or PS-AH34-1).
- Since the Hybrid Base™ Bridge and Prosthetic Screw are small they must be handled with caution to avoid the patient from swallowing or inhaling them.
- The Hybrid Base™ Bridge and Prosthetic Screw have not been evaluated for safety, heating, migration or compatibility in the Magnetic Resonance Imaging (MRI) environment.
- Place implant-borne restorations only in occlusion when the implant is completely osseointegrated.
- Always place temporary restorations out of occlusion.
- Allergies to the titanium alloy grade 5 ELI (TiAl<sub>6</sub>V<sub>4</sub> ELI), or contents of the alloy may very rarely occur.
- The Elos Accurate® API® Bridge components must only be used and handled by dental professionals.
- The use of torque value higher than the recommended for the implant system according to the manufacturer's instruction may result in damage to the Hybrid Base™ Bridge, the Prosthetic screw and/or the implant. The use of torque values lower than those recommended for the implant system according to the manufacturer's instruction may result in loosening of the Hybrid Base™ Bridge and/or the Prosthetic screw.
- For any lab work of the Hybrid Base™ Bridge, the included Process screw must be used.
- The included Prosthetic screw may only be used to permanently fix the finished restoration.

#### 6 Compatibility information

The Elos Accurate® API® Bridge product is available for a variety of implant platforms and sizes. The individual product label refers to which implant platform and size the product is compatible to.

#### 7 Cleaning and disinfection

The Elos Accurate® API® Bridge prosthetic components are delivered non-sterile. Prior to attachment of the prosthetic restoration in the patient's mouth, it must be cleaned and disinfected followed by sterilization. Following cleaning and disinfection method is recommended:

1. Clean rinsing under flowing water while brushing outer and inner side with adequate brushes.
2. The pre-treated restoration can be cleaned by use of compatible cleaning solution (e.g. Cidezyme/Enzol, Johnson & Johnson).
3. The cleaned restoration can be disinfected in compatible disinfectant (e.g. Cidex OPA, Johnson & Johnson).
4. If the cleaned and disinfected prosthetic restoration has to be packed, it should be packed in adequate packaging suitable for steam sterilization.

Cleaning and disinfection should be processed as specified by the manufacture. Cleaned and disinfected parts should be packed in adequate packaging according to instruction for use supplied by the bridge material manufacturer.

#### 8 Sterilization

The Elos Accurate® API® Bridge components are delivered non-sterile. The recommended sterilization procedure is full cycle pre-vacuum steam sterilization at a temperature of 132 °C (270° F) for 4 minutes, 4 pulse minimum or full cycle gravity steam sterilization at 132°C (270°F) for 15 minutes minimum. The products are ready for use 20 minutes after autoclaving process.

**Caution:** The device must be used immediately after sterilization. Do not store sterilized devices.

#### 9 Procedure

##### 9.1 Use and handling by dental technician

Gently open the outer packaging and take out the blister box. By use of a scissor, separate the Prosthetic screw from the Hybrid Base™ Bridge, Process screw and Model analog. Make sure to separate the blister box in the dotted line in order to protect the sealing of the Prosthetic screw. The Prosthetic screw must be provided to the dental clinic with the prosthetic restoration.

##### 9.1.1 Manufacture the prosthetic restoration

**Recommended bridge material:** 3M LAVA™ Frame zirconia

**Compatible design software:** 3Shape, Exocad, and Dental Wings

##### Caution:

- Ensure to follow the selected dental bridge restoration material minimum wall thickness and parameters according to the ceramic material supplier's instruction for use.
- The dental bridge must be manufactured with milling equipment adequate for the designed geometry
- A precise milling of the geometry requires millers of 1 mm in diameter or smaller.
- The dental bridge should be sintered in accordance with the instruction issued by the ceramic material supplier.

##### Using digital workflow (intra oral scanning):

1. Import the digitized patient situation from the dental clinical intraoral scan to the design software by use of Elos Accurate® Scan Body in order to identify the position and orientation of the respective implant.
2. Import library file from Elos Medtech webpage and select relevant implant platform from the library in order to facilitate precise design of

inner geometry of the dental bridge interfacing the Hybrid Base™ Bridges.

3. Design dental bridge in the design software or scanned teeth set up validated from the clinic in the patient's mouth.
4. Export STL file to milling equipment. The STL file contains the required milling template for the inner geometry of the prosthetic restoration.
5. Process and finalize the prosthetic restoration in milling equipment according to the manufacturer's instruction for use.
6. Create a digital working model in the design software.
7. Export the STL file from the design software and send the STL file to your 3D printer or external 3D print provider.
8. Place included Model analogs in 3D printed working model.

#### Using semi-digital workflow (desktop scanning):

1. Create working model and place the enclosed model analog.
2. Place preferable Elos Accurate Scan Body in the model analog.
3. By use of an Elos Accurate Scan Body, scan the working model.
4. Remaining steps as for the digital workflow.

#### 9.1.2 Preparing the prosthetic restoration for cementing

1. Blast inner geometry of the milled dental bridge interfacing the Hybrid Base™ Bridge with aluminum oxide  $\leq 50 \mu\text{m}$  and blasting pressure of 2 bar.
2. Before cementing, clean the surface thoroughly with alcohol.

Note: Blasting of Hybrid Base™ Bridges is not necessary.

#### 9.1.3 Cementing the prosthetic restoration

1. Place Hybrid Base™ Bridge in the milled prosthetic restoration with correct alignment of screw channels. The Hybrid Base™ Bridges are kept in place by the GuideGrip™ Technology (retention fingers).
2. Check the clearance of the screw channels by inserting the process screws in the milled dental bridge / Hybrid Base™ Bridges.
3. Place dental bridge assembly on the working model with Model analogs and check fit. Hand-tighten Process screws. Always use Elos Accurate® Prosthetic Screwdriver for hexalobular screws supplied by Elos Medtech.
4. Remove the milled dental bridge from the working model (the Hybrid Base™ Bridges are kept in the working model by the process screws).
5. Seal the screw channels with wax.
6. Apply self-adhesive dental cement on the the Hybrid Base™ Bridges. Use only suitable self-adhesive cementation systems for the dental bridge material used. Follow the instruction for use of both the dental bridge material and cement material manufacturer.

**Recommended cement:** 3M ESPE RelyX™ self-adhesive cement

1. Bond the milled dental bridge to the Hybrid Base™ Bridges placed in the working model. The milled dental bridge must fit the Hybrid Base™ Bridges without gaps and voids.
2. Immediately remove excess cement from the Hybrid Base™ Bridges.

**Caution:** Include the Prosthetic screws in its unbroken packaging when sending the final dental bridge restoration to the dental clinic.

#### 10 Use and handling for the dentist

The dentist receives the final dental bridge restoration / working model with the Prosthetic screws in unbroken packaging from the dental lab.

1. Identify and unpack the Prosthetic screws.
2. Clean, disinfect and sterilize prosthetic restoration and Prosthetic screws according to section 7 and 8 in this instruction for use.
3. Remove the healing cap, closure screw or temporary restoration from patient's mouth.
4. Gently, insert the bridge restoration into the patient's mouth in proper position to the implants.
5. Place corresponding Prosthetic screws in the bridge restoration and tighten screws with Elos Accurate® Prosthetic Screwdriver for

hexalobular screws, using torque values specified by the implant manufacturer.

6. Close screw channels before cementing (e.g. with wax). This allows a later removal of the Hybrid Base™ Bridge prosthetic components in case a bridge replacement is necessary.
7. Additionally, if cementing of the custom bridge and the Hybrid Base™ Bridges has to be done intraorally, place Hybrid Base™ Bridges in the milled prosthetic restoration with correct rotation of screw channels.
8. Seal the screw channels with wax.
9. Apply self-adhesive dental cement on the the Hybrid Base™ Bridges. Use only suitable self-adhesive cementation systems for the dental bridge material used. Follow the instruction for use of both the dental bridge material and cement material manufacturer.

**Recommended cement:** 3M ESPE RelyX™ self-adhesive cement

10. Bond the milled dental bridge restoration to the Hybrid Base™ Bridges. The milled dental bridge must fit the Hybrid Base™ Bridges without gaps and voids. Immediately remove excess cement from the Hybrid Base™ Bridges.

#### 11 Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

#### 12 Validity

Upon publication of this instruction for use, all previous versions are superseded.

#### 13 Storage and Handling:

Elos Accurate® API® Bridge should be stored at room temperature.

#### 14 Disposition:

The prosthetic restoration must be disposed as biological waste.

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