Elos Prosthetic Screw

Instruction for use - Инструкция за
употреба - Návod k použití - Brugsanvisning
- Gebrauchsanleitung - Kasutusjuhend -
Οδηγίες χρήσης - Instrucciones de uso -
Mode d'emploi - Upute za upotrebu -
Notkunarleiðbeiningar - Istruzioni per l'uso -
Lietošanas norādījumi - Naudojimo
instrukcija - Használati útmutató -
Gebruiksaanwijzing - Bruksanvisning -
Instrukcja użytkowania - Instruções de uso -
Instruțiuni de utilizare - Návod na použitie -
Navodila za uporabo - Käyttöohjeet -
Bruksanvisning


List of content

English – Instruction for use.....	3
Български – Инструкция за употреба (Bulgarian)	4
Čeština – Návod k použití (Czech).....	5
Dansk – Brugsanvisning (Danish).....	6
Deutsch – Gebrauchsanleitung (German).....	7
Eesti keel – Kasutusjuhend (Estonian).....	8
Ελληνικά — Οδηγίες χρήσης (Greek)	9
Español – Instrucciones de uso (Spanish)	10
Français – Mode d’emploi (French)	11
Hrvatski – Upute za upotrebu (Croatian).....	12
Íslenska – Notkunarleiðbeiningar (Icelandic).....	13
Italiano – Istruzioni per l’uso (Italian).....	14
Latviski — Lietošanas norādījumi (Latvian).....	15
Lietuvių k. – Naudojimo instrukcija (Lithuanian)	16
Magyar – Használati útmutató (Hungarian).....	17
Nederlands – Gebruiksaanwijzing (Dutch).....	18
Norsk – Bruksanvisning (Norwegian).....	19
Polski – Instrukcja użytkowania (Polish)	20
Português – Instruções de uso (Portuguese).....	21
Română – Instrucțiuni de utilizare (Romanian)	22
Slovenčina – Návod na použitie (Slovenian).....	23
Slovenščina – Navodila za uporabo (Slovak).....	24
Suomi – Käyttöohjeet (Finnish).....	25
Svenska – Bruksanvisning (Swedish)	26

English – Instruction for use

1. **Indications for Use:** The Elos Prosthetic Screw is intended for fixation of prosthetic restorations to endosseous dental implants
2. **Device Description:** The Elos Prosthetic Screw is manufactured from biocompatible titanium alloy grade 4/5 ELI (TiAl6V4 ELI). Depending on the specific variant, the Elos Prosthetic screw includes biocompatible coatings.
3. **Contraindications:** Allergies to the alloy, or contents of the alloy may very rarely occur.
4. **Operating procedure:** Ensure that the prosthetic restoration is positioned correctly on the implant. Tighten the Elos prosthetic screw to the recommended torque using an appropriate tool. Recommended torque value is stated on label.
5. **Warnings:**
 - The use of torque values higher than the recommended may result in damage to the prosthetic restoration, the screw and/or the implant. The use of torque values lower than those recommended may result in loosening of the prosthetic restoration.
 - The Elos Prosthetic Screw is for single use only. Reuse can result in loss of functionality and/ or infections.
 - Since the Elos Prosthetic Screw is small it must be handled with caution to avoid the patient swallowing or inhaling them.
6. **Precautions:** The Elos Prosthetic Screw must only be used and handled by dental professionals who have approved training in dentistry.
For any lab work and any try-in of the prosthetic restoration, a lab screw must be used. The Elos Prosthetic Screw may only be used to permanently fix the finished restoration.
7. **Cleaning and Sterilization:** Prior to installing, the Elos Prosthetic Screw must be cleaned with water and brush followed by disinfection and sterilization. The recommended sterilisation procedure is a full cycle pre-vacuum steam sterilization at a temperature of 134 °C (273 °F) for 3 minutes minimum or 121°C (250°F) for 15 minutes minimum.
8. **Compatibility:** For compatibility of the Elos Prosthetic Screw please refer to information on individual product label.
9. **MRI Statement:** The Elos Prosthetic Screw has not been evaluated for safety, heating, migration or compatibility in the Magnetic Resonance Imaging (MRI) environment.
10. **Storage and Handling:** The Elos Prosthetic Screw should be stored at room temperature protected from direct sunlight.
11. **Disposal:** Must be disposed as biological waste.
12. **Symbols:**

	Catalogue number		Prescription only
	Batch code		Do not use if package is damaged
	Manufacturer		Non sterile
	Consult instructions for use		Recommended torque
	Do not re-use		

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