Elos Torque Wrench & Driver

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
Elos Torque Wrench & Driver

1 Intended use

The Elos Torque Wrench & Driver is adaptable to a variety of tools in order to apply torque for either tightening/loosening/removing a screw and or implant during a surgical procedure.

2 Product Description

The Elos Torque Wrench (main body, head and Elos Driver) is made of lightweight titanium and stainless steel. The Elos Torque Wrench & Driver product design, finish and quality make the tool well suited for placement and adjustment of implants, abutments and screws during surgery and other clinical procedures. The Elos Torque Wrench & Driver is easy to prepare for cleaning, disinfection and sterilization.

Exchangeable screwdrivers and the Elos Driver make the tool compatible to various implant systems.

3 Use

The Elos Torque Wrench & Driver is delivered non-sterile. Prior to use it must be cleaned with water and brush followed by sterilization according to section 5 and 6. The use of the Elos Torque Wrench & Driver is described below:

1. After cleaning make sure the first line on the scale of the Elos Torque Wrench arm aligns with the arrow (*see figure 1.A*).
2. Attach the head to the body of the Elos Torque Wrench by pushing the components together and turning them in opposite direction until an audible click is achieved (*see figure 2*).
3. The arrow on the head of the Elos Torque Wrench shows the direction in which the Elos Torque Wrench & Driver is functioning, and the head is labeled "IN" (*see figure 3.A*) and "OUT" (*see figure 3.B*).
"IN" denotes clockwise rotation; "OUT" denotes counterclockwise rotation.
In order to change direction flip the Elos Torque Wrench to the opposite side (*see figure 3*) and push the Elos Driver into the head of the Elos Torque Wrench until an audible click is achieved (*see figure 4*).
4. Insert the tool into the Elos Driver.
5. Place the tool in the implant/screw and turn the Elos Torque Wrench arm in the direction of the arrow until the desired torque is archived (*see figure 6*). The applied torque must be in accordance with the screw/implant manufacturers' instructions.
6. After each use clean and sterilize the Elos Torque Wrench & Driver according to section 5 and 6.

4 Warnings and precautions

Potential adverse occurrences may include but are not limited to the following list. Proper placement and use of this instrument may help to mitigate these risks.

- **If overloaded or dropped**

If overloaded, dropped or in other way mishandled, the product must no longer be used since correct functionality can no longer be guaranteed.

- **Do not use hydrogen peroxide**

This product must be sterilized prior to use and it must not be cleaned by using hydrogen peroxide.

- **Scale alignment**

The Elos Torque Wrench & Driver does not require recalibration but before each use, make sure that the Elos Torque Wrench & Driver is fully functional and that the first line on the scale aligns with the arrow (*see figure 1.A*).

- **Ensure correct assembly**

Before each use, make sure that all parts are correct assembled (no loose parts are allowed).

- **Do not exceed max scale**

The arm of the Elos Torque Wrench must not go beyond the end of the scale (*see figure 6.B*), as this could result in permanent deformation of the Elos Torque Wrench arm i.e. future inaccurate torque readings (*see figure 1.B*).

- **No scale indication when used as ratchet wrench**

If the Elos Torque Wrench & Driver is used as a ratchet wrench it does NOT provide the user any indication of what torque that is used.

Notice: the torque applied must not exceed 150Ncm when the Elos Torque Wrench & Driver is used as a ratchet wrench.

- **Use screw/screwdriver/Driver/implant manufacturer torque recommendations**

Please check the individual screw/screwdriver/Elos Driver/implant manufacturers' instruction for use or label regarding the maximum allowed torque. Do not under any circumstances apply higher torque than stated even if the Elos Torque Wrench & Driver is used as a ratchet wrench (without using the Elos Torque Wrench arm).

- **Cleaning – do not use metal brush or steel wool**

Do not use a metal brush or steel wool to clean the Elos Torque Wrench or Driver as this will damage the instrument.

- **Only for use by professionals**

The Elos Torque Wrench & Driver must only be used by professionals.

5 Cleaning and disinfection

The following cleaning and disinfection method is recommended:

1. The body and head of the Elos Torque Wrench are assembled by delivery. To remove the head from the body of the Elos Torque Wrench press into the recess area of the head/main body and gently pull the two parts from each other.
If the Elos Driver is attached to the head of the Elos Torque Wrench the Elos Driver must be removed before the head and body of the Elos Torque Wrench are disassembled.
Press the Elos Driver to remove it from the head of the Elos Torque Wrench (*see figure 5*).
2. The three separated parts are now ready to be cleaned with water and brush and sterilized according to section 6.

6 Sterilization

The Elos Torque Wrench & Driver is delivered non-sterile. Steam sterilize with a pressure of 2 atm. for 15 min. at 121°C or with a pressure of 3 atm. for 3 min. at 134°C, according to EN ISO 17665-1:2006 is recommended.

Caution: Accompanying packaging material CANNOT be autoclaved unless it is explicitly stated that the material is suitable for sterilization.

7 Further information

For additional information about the use of the Elos Torque Wrench & Driver, please contact your local sales representative.

8 Validity

Upon publication of this instruction for use, all previous versions are superseded.

9 Storage and Handling

The Elos Torque Wrench & Driver should be stored at room temperature. Store at 10°- 40°C (50°-104°F).

10 Disposal

The Elos Torque Wrench & Driver must be disposed as biological waste.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

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