# Elos Drill Extender

Instruction for use - Upute za upotrebu -  
Brugsanvisning - Gebruiksaanwijzing - Käyttöohjeet -  
Mode d'emploi - Gebrauchsanleitung - Használati  
útmutató - Istruzioni per l'uso - Bruksanvisning -  
Instruções de uso - Instrucțiuni de utilizare -  
Instrucciones de uso - Anvisningar för användningen

## List of content

English – Instruction for use .....	4
Hrvatski/ Croatian – Upute za upotrebu .....	5
Dansk/ Danish – Brugsanvisning .....	6
Nederlands/ Dutch – Gebruiksaanwijzing .....	7
Français/ French – Mode d’emploi .....	9
Deutsch/ German – Gebrauchsanleitung .....	10
Magyar/ Hungarian – Használati útmutató .....	11
Italiano/ Italian – Istruzioni per l’uso .....	12
Norsk/ Norwegian – Bruksanvisning .....	13
Português/ Portuguese – Instruções de uso .....	14
Română/ Romanian – Instrucțiuni de utilizare .....	15
Srpski/ Serbian – Uputstva za upotrebu .....	16
Español/ Spanish - Instrucciones de uso .....	17
Svenska/ Swedish – Bruksanvisning .....	18

If package is Broken - DO NOT USE the product - Nemojte upotrebljavati ako je pakiranje oštećeno – Anvend ikke hvis emballagen er beskadiget - Niet gebruiken indien de verpakking beschadigd is - Älä käyttää, jos pakkaus on vahingoittunut - Ne l'utilisez pas si l'emballage est endommagé - Nicht verwenden, wenn die Verpackung beschädigt ist - Ne használja, ha a csomagolás sérült - Non utilizzare se la confezione è danneggiata - Produktet skal ikke brukes hvis emballasjen er skadet - Não utilizar se a embalagem estiver danificada - A nu se utiliza dacă ambalajul este deteriorat - No utilizar si el paquete está dañado - Använd inte om förpackningen är skadad

Figure 1.



Figure 2

## English – Instruction for use

### Elos Drill Extender

#### 1 Intended use

The intended use of the Elos Drill Extender is to extend the distance between a dental drill or screwdriver and the dental hand piece or Elos Torque Wrench during surgery and prosthetic procedures in the oral cavity.

#### 2 Product Description

The Elos Drill Extender is for extending tools with ISO 1797 connections. The device is made of stainless steel. The Elos Drill Extender design, finish and quality makes the tool well suited for extending a drill or dental tool and to placement and adjustment of implants, abutments and prosthetic screws. The Elos Drill Extender is easy to prepare for cleaning, disinfection and sterilization.

#### 3 Use

The Elos Drill Extender is delivered non-sterile. Prior to use the device must be cleaned with water and brush followed by sterilization according to section 5 and 6. Use the Elos Drill Extender as described below:

1. Figure 1: After cleaning attach the Elos Drill Extender (A) to the Elos Driver ISO connection (B) or into the dental hand piece ISO connection (C). Push the Elos Drill Extender into the ISO connection until firmly seated.
2. Figure 2: Attach the desired dental tool or drill (D) into the Elos Drill Extender (A), and make sure that the tool is firmly seated in the Elos Drill Extenders ISO connection.

#### 4 Warnings and precautions

Potential adverse occurrences may include but are not limited to the following list. Proper placement and use of this instrument may help to mitigate these risks.

- **If overloaded**

If overloaded, dropped or in other way mishandled, the device must no longer be used since correct functionality can no longer be guaranteed.

- **Damaged device**

If the Elos Drill Extender is damaged it should be disposed, as described in section 10.

- **Do not use hydrogen peroxide**

This device must not be cleaned using hydrogen peroxide.

- **Ensure correct assembly**

Before use, make sure that all parts are correctly assembled.

- **Must be stored unattached**

The Elos Drill Extender must be stored unattached to other instruments.

- **Max torque 100Ncm**

Never apply torque > 100Ncm to the Elos Drill Extender, as this may damage the ISO connection. **However never exceed the screw/implant manufacturer torque recommendations**

- **Max speed 3000 rev/min**

Never use drilling speed > 3000 rev/min, as this can damage the Elos Drill Extender.

- **Cleaning**

The Elos Drill Extender must always be cleaned using brush and water.

- **Only for use by professionals**

The device must only be used by professionals.

#### 5 Cleaning and disinfection

The following cleaning and disinfection method is recommended:

1. Remove the tool from the Elos Drill Extender by gently pulling the two parts from each other.
2. Remove the Elos Drill Extender from Elos Torque Wrench or drilling device by gently pulling the two parts from each other.
3. The Elos Drill Extender is now ready to be cleaned with water and brush and sterilized according to section 6.

#### 6 Sterilization

The Elos Drill Extender is delivered non-sterile. Steam sterilize with a pressure of 2 atm. for 15 min. at 121°C or with a pressure of 3 atm. for 3 min. at 134°C, according to EN ISO 17665-1:2006 is recommended.

#### 7 Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

#### 8 Validity

Upon publication of this instruction for use, all previous versions are superseded.

#### 9 Disposal

Device must be disposed as biological waste.

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