Elos Accurate® Scan Body

Instruction for use - Инструкция за
употреба – Návod k použití –
Brugsanvisning – Gebrauchsanleitung –
Kasutusjuhend – Instrucciones de uso –
Οδηγίες χρήσης – Mode d'emploi – Upute
za upotrebu – Istruzioni per l'uso –
Lietošanas norādījumi – Naudojimo
instrukcija – Használati útmutató –
Gebruiksaanwijzing – Bruksanvisning –
Instrukcja użytkowania – Instruções de uso
– Instrucțiuni de utilizare – Návod na
použitie – Navodila za uporabo –
Käyttöohjeet – Bruksanvisning – Kullanma
Talimatı

Disclaimer of liability: Elos Medtech Pinol products may only be used according to the instructions and recommendations. The user of Elos Medtech Pinols products should determine their suitability for particular patients and indications. Elos Medtech Pinol disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or damages arising out of or in connection with any errors in professional judgement or practice in the use or placement of the Elos Medtech Pinol products.

Please note that not all products are regulatory cleared/released for sale in all markets even though the language is present in this instruction.

List of content

English - Instruction for use (English)	5
Български – Инструкция за употреба (Bulgarian)	7
Hrvatski – Upute za upotrebu (Croatian).....	9
Čeština - Návod k použití (Czech)	11
Dansk – Brugsanvisning (Danish).....	13
Nederlands – Gebruiksaanwijzing (Dutch).....	15
Eesti keel – Kasutusjuhend (Estonian).....	17
Suomi – Käyttöohjeet (Finnish).....	19
Français - Mode d'emploi (French)	21
Deutsch - Gebrauchsanleitung (German)	23
Ελληνικά - Οδηγίες χρήσης (Greek).....	25
Magyar - Használati útmutató (Hungarian)	27
Italiano - Istruzioni per l'uso (Italian)	29
Latviski - Lietošanas norādījumi (Latvian).....	31
Lietuvių k. - Naudojimo instrukcija (Lithuanian)	33
Norsk – Bruksanvisning (Norwegian).....	35
Polski - Instrukcja użytkowania (Polish)	37
Português - Instruções de uso (Portuguese).....	39
Română - Instrucțiuni de utilizare (Romanian)	41
Slovenčina - Návod na použitie (Slovak).....	43
Slovenščina – Navodila za uporabo (Slovenian)	45
Español - Instrucciones de uso (Spanish).....	47
Svenska – Bruksanvisning (Swedish)	49
Türkçe - Kullanma Talimatı (Turkish).....	51

**LOT**

Batch code - Партиден номер - Šifra serije - Kód šarže - Batchkode - Batchcode - Partii kood - Eräkoodi - Code du lot - Chargencode - Κωδικός παρτίδας - Tételszám - Numero di lotto - Partijas kods - Partijos kodas - Batchkode - Numer partii - Código de lote - Cod lot - Číslo šarže - Koda serije - Código de lote - Batchkod - Parti kodu -

REF

Catalogue number - Каталоген номер - Kataloški broj - Katalogové číslo - Katalognummer - Catalogusnummer - Katalooginumber - Luettelonumero - Référence catalogue - Katalognummer - Αριθμός καταλόγου - Katalógusszám - Codice catalogo - Kataloga numurs - Katalogo numeris - Katalognummer - Numer katalogowy - Número do catálogo - Număr catalog - Katalógové číslo - Kataloška številka - Número de catálogo - Katalognummer - Katalog numarası

Manufacturer - Производител - Proizvođač - Výrobce - Producent - Fabrikant - Tootja - Valmistaja - Fabricant - Hersteller - Κατασκευαστής - Gyártó - Produttore - Izgatavotājs - Gamintojas - Produsent - Producent - Fabricante - Producător - Výrobca - Proizvajalec - Fabricante - Tillverkare - Üretici

Date of manufacture - Дата на производство - Datum proizvodnje - Datum výroby - Fremstillingsdato - Productiedatum - Tootmiskuupäev - Valmistuspäivämäärä - Date de fabrication - Productiedatum - Ημερομηνία κατασκευής - Gyártás dátuma - Data di produzione - Pagaminimo data - Izgatavošanas datums - Produksjonsdato - Data produkcji - Data de fabrico - Data fabricației - Dátum výroby - Data fabricației - Fecha de fabricación - Tillverkningsdatum - Üretim tarihi

By prescription only - Само по предписание - Isključivo na recept - Pouze na předpis - Receptpligtig - Alleen op recept - Ainult meditsiiniliseks kasutuseks - Vain lääkärin määräyksestä - Sur ordonnance uniquement - Verschreibungspflichtig - Μόνο κατόπιν εντολής ιατρού - Kizárolag receptre kapható - Solo su prescrizione - Tikai pēc norādījuma - Receptinis - Kun etter forordning - Wyłącznie z przepisu lekarza - Apenas com receita médica - A se elibera numai pe bază de prescripție medicală - Len na predpis - Samo na recept - Solamente con receta - Endast med recept - Yalnızca reçeteyle satılır

Consult instructions for use - Направете справка с инструкциите за употреба - Pročítajte upute za uporabu - Prostudujte návod k použití - Se brugsanvisning før brug - Raadpleeg gebruiksaanwijzing - Lugege kasutusjuhendit - Lue käyttöohjeet ennen käyttöä - Consulter le mode d'emploi - Gebrauchsanleitung beachten - Συμβουλευτείτε τις οδηγίες χρήσης - Olvassa el a használati útmutatót - Consultare le istruzioni per l'uso - Skatiet lietošanas instrukcijas - Žr. naudojimo instrukcijas - Les bruksanvisningen - Zapoznać się z instrukcją użytkowania - Consulte as instruções de uso - Consultați instrucțiunile de utilizare - Prečítajte si návod na použitie - Glejte navodila za uporabo - Consulte las instrucciones de uso - Läs bruksanvisningen - Kullanna talimatlarini inceleyin -

Do not use if package is damaged - Не използвайте, ако целостта на опаковката е нарушена - Nemojte upotrebljavati ako je pakiranje oštećeno - Nepoužívejte, je-li obal poškozený - Må ikke anvendes, hvis pakken er beskadet - Niet gebruiken indien de verpakking beschadigd is - Ärge kasutage kahjustatud pakendit - Älä käytä, jos pakkaus on vahingoittunut - Ne l'utilisez pas si l'emballage est endommagé - Nicht verwenden, wenn die Verpackung beschädigt ist - Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά - Ne használja, ha a csomagolás sérült - Non utilizzare se la confezione è danneggiata - Nelietot iepakojumu, ja tas ir bojāts - Nenaudokite, jei pakuotė pažeista - Produktet skal ikke brukes hvis emballasjen er skadet - Nie używać w razie uszkodzenia opakowania - Não utilizar se a embalagem estiver danificada - A nu se utiliza dacă ambalajul este deteriorat - Nepoužívajte, ak je poškodený obal - Ne uporabljajte, če je embalaža poškodovana - No utilizar si el paquete está dañado - Använd inte om förpackningen är skadad - Hasar görmüş paketleri kullanmayın

Medical device - Медицинско изделие - Medicinsko sredstvo - Lékařské zařízení - Medicinsk udstyr - Medisch apparaat – Meditsiiniseade - Lääketieteellinen laite - Dispositif médical - Medisch apparaat - Ιατροτεχνολογικό προϊόν - Orvosi eszköz - Dispositivo medico - Medicinos priemonė - Medicinas ierīce - Medisinsk utstyr - Wyrób medyczny - Dispositivo médico - Dispozitiv medical - Zdravotnícka pomôcka - Dispozitiv medical - Dispositivo médico - Medicinsk utrustning - Tibbi Cihaz

Unique Device Identifier - Унікален ідентифікатор на изделието - Jedinstveni identifikator uređaja - Jedinečný identifikátor zdravotníckého prostriedku - Unik anordningsidentifikator - Unieke apparaat-ID - Kordumatu identifitseerimistunnus - Yksilöllinen laitetunnus - Identifiant unique de dispositives - Eindeutige Produktkennung - Μοναδικό αναγνωριστικό συσκευής - Egyedi eszközazonosító - Identificatore di dispositivi univoci - Ierīces unikālais identifikators - Unikalus prietaiso identifikatorius - Unik enhetsidentifikator - Unikatowy identyfikator urządzenia - Identificador único de dispositivo - Identificatorul unic al dispozitivului - Jedinstveni identifikator uređaja - Jedinečný identifikátor pomôcky - Enolični identifikator pripomočka - Identificador de producto único - Unik produktidentifizierung - Benzersiz Cihaz Tanımlayıcı

Non sterile – Нестерилен - Nije sterilno – Nesterilni - Ikke-steril – Niet-steriel - Pole steriilne – Epästeriili - Non stérile – Unsteril - Μη αποστειρωμένο - Nem steril - Non sterile – Nesterilis – Nesterilus - Ikke-sterilt – Niejałowy - Não esterilizado - Prodnis nesterilizat – Nesterilné - Ni sterilno - No estéril - Ikke-steril - Non-steril

Keep dry – Pazete Пазете сухо – Držite na suhom – Udržujte v suchu – Hold tørt – Droog bewaren – Hoida kuivana – Säilytä kuivana – Tenir au sec – Vor Feuchtigkeit schützen – Διατηρείτε στεγνό– Tartsa szárazon – Conservare in luogo asciutto – Turēt sausā vietā – Laikyti sausai– Holdes tørr – Przechowywać w suchym miejscu – Mantenha seco – A se menține uscat – Uchovávaťe v suchu – Hranite na suhem – Mantener en un lugar seco – Förvara torrt – Kuru yerde saklayın

Keep away from sunlight - Пазете далече от слънчева светлина - Držite podalje od Sunčeve svjetlosti - Chraňte před slunečním světlem - Hold ude af sollys - Uit de buurt van zonlicht houden - Kaitsta päikesevalguse eest - Säilytä suojattuna auringonvalolta - Tenir à l'abri de la lumière du soleil - Vor Sonnenlicht schützen - Διατηρείτε μακριά από το ηλιακό φως - Napfénytől távol tartandó - Conservare al riparo dalla luce del sole - Sargāt no saules gaismas - Saugoti nuo saulės šviesos - Holdes vekk fra sollys - Chronić przed światłem słonecznym - Mantenha afastado da luz solar - A se feri de lumina soarelui - Uchovávaťe mimo slnečného svetla - Ne izpostavljajte sončni svetlobi - Mantener alejado de la luz solar - Háll borta från solljus - Güneş ışığından uzak tutun

English - Instruction for use (English)

Elos Accurate® Scan Body

1 Intended purpose

The Elos Accurate® Scan Body is intended to be used together with a dental scanner to locate the exact position of an implant/analog or abutment in relation to prosthetic restorations. The Scan Body is intended to be mounted directly on the implant/analog or abutment using the Elos Accurate® Scan Body Driver. The Scan Body is intended to be used directly in a patient's mouth when using an intra oral scanner, or on a model of the patients' teeth setup when using a desktop scanner.

2 Product Description

The Scan Body is manufactured from biocompatible titanium and PEEK dependent on the variant. The corresponding Scan Body Driver used to tighten the integrated screw is manufactured from biocompatible stainless steel. The product is available for a variety of implant platforms. For specific product descriptions, please refer to individual product labels.

3 Chemical composition

Chemical composition of materials used for the Scan Body and Scan Body Driver is presented in table 1-4.

Table 1: Chemical composition of the ScanBody22™ made of AccuScan® coated titanium, titanium alloy and anodized titanium alloy.

Material	Composition	Patient contacting	CAS nr.
Anodized Titanium alloy (Color oxidized)	33.3 at.% Ti + 66.6 at.% O (TiO ₂)	Yes (Outer layer)	1317-80-2
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	Yes (Pure Metal)	99906-66-8
AccuScan®	33.3 at.% Ti + 66.6 at.% O (TiO ₂)	Yes (Outer layer)	1317-80-2
Titanium Grade 4	100 wt.% Ti	No (Inner layer)	7440-32-6

Table 2: Chemical composition of the Scan Body made of Titanium alloy, anodized titanium alloy and PEEK.

Material	Composition	Patient contacting	CAS nr.
Anodized Titanium alloy (Color oxidized)	33.3 at.% Ti + 66.6 at.% O (TiO ₂)	Yes (Outer layer)	1317-80-2
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	Yes (Pure Metal)	99906-66-8
PEEK (polyetheretherketone)	C ₁₉ H ₁₄ O ₃	Yes	29658-26-2

Table 3: Chemical composition of the Scan Body Driver made of Stainless steel.

Material	Composition	Patient contacting	CAS nr.
Sandvik 4C27A	84.8 wt.% Fe, 13 wt.% Cr, 1.2 wt.% Mo, 0.8 wt.% Ni, 0.2 wt.% S	Yes (Pure Metal)	65997-19-5
AISI 304 L	67-71.5 wt.% Fe, 17.5-19.5 wt.% Cr, 8.0-10.5 wt.% Ni, 2 wt.% Mn, 1 wt.% Si	Yes (Pure Metal)	65997-19-5

Table 4: Chemical composition of the ScanBody22™ Driver made of Stainless steel.

Material	Composition	Patient contacting	CAS nr.
Sandvik 1RK91	71,7 wt.% Fe, 12 wt.% Cr, 9 wt.% Ni, 4 wt.% Mo, 2 wt.% Cu, 0.9 wt.% Ti and 0.4 wt.% Al	Yes (Pure Metal)	65997-19-5

4 Basic UDI

Scan Body: 5712821004001UP
 Scan Body Driver: 5712821004003UT
 Desktop Scan Body: 5712821004002UR

5 Indication for use

The Scan Body is intended to be mounted directly on the implant/analog or abutment using the Scan Body Driver. The Scan Body is intended to be used directly in a patient's mouth when using an intra oral scanner, or on a model of the patient's teeth setup when using a desktop scanner.

6 Intended patient population

The Scan Body is intended to be used on patients that are healthy enough to have a dental implant inserted in the jaw.

7 Intended users

The Scan Body should only be used by dentists or dental technicians who are experienced within dental implantology.

Use environment

Dental laboratories and dental clinics.

8 Clinical benefits

The Scan Body does not have a clinical benefit. It is simply used to successfully assist as a scanning tool to support an optimal dental restoration.

9 Expected lifetime

Expected lifetime of up to 1 year or 100 use cycles.

10 Contraindications

None known

11 Warnings and precautions

- The Scan Body may only be used with the original branded implant/analog or abutment stated on the product label.
- Make sure that the seating surface, upon which the scan body is to be placed, is clean and free of residue and liquids as well as bone or soft tissue when mounting.
- Make sure that the Scan Body is not in contact with surrounding teeth when mounted on the implant, analog or abutment.
- Since the Scan Body is small it must be handled with caution to avoid the patient from swallowing or inhaling it.
- Tighten the integrated screw must only be done slightly by hand using the Scan Body Driver (max. 5Ncm). To high tightening torque can damage the device which may cause inaccurate scan results.
- DO NOT use a Scan Body that has been damaged or modified in any way. Use of a damaged or modified Scan Body may cause inaccurate scan result.
- Do not try to disassemble the Scan Body.
- Do not try to rotate the Scan Body, when it is attached to the implant.
- The Scan Body may only be used with a compatible digital library file installed in the design software.
- The Scan Body must only be used and handled by dental professionals.
- Allergies to the materials or contents of the materials may very rarely occur.
- Inspect Scan Body/Scan Body Driver under normal lighting for the removal of visible soil and for damage or wear: Any Scan Body/Scan Body Driver showing visible damage or wear must be disposed of according to section 18.

12 Compatibility information

The Scan Body is available for a variety of implant platforms and sizes. The individual product label refers to which implant platform and size the product is compatible with.

13 Cleaning and Sterilization

The Scan Body is delivered non-sterile. Prior to attachment of the Scan Body in the patient's mouth, it must be cleaned followed by sterilization. The Scan Body shall be cleaned and sterilized before use (initially as well as between use). Further instruction can be found in "Cleaning and sterilization guideline" at <https://elosmedtech.com/IFU/>.

- Clean according to the guidance in the "Cleaning and sterilization guideline".

2. Sterilize according to the guidance in the "Cleaning and sterilization guideline". Summarized below.

Note: For US: Steam Sterilization Cycle 132°C (270°F)/4min.

Procedure	Dynamic-air-removal sterilization Cycle	
Exposure time	4 min.	3 min.
Temperature	132°C (270°F)	134° C (273°F)
Drying time	20 min.	

Caution:

- Do not clean the Scan Body with ultrasound as this can damage the product.
- As the Scan Body is a high precision tool, it may not be used while still warm from the autoclave, as the thermal expansion may lead to inaccurate scanning. Before using the Scan Body, it must be cooled to 18-37°C in temperature.
- Accompanying packaging material CANNOT be autoclaved unless it is explicitly stated that the material is suitable for sterilization.

14 Procedure

14.1 Scanning

1. Clean and sterilize the Scan Body according to the procedure described in section 13. Sterilization is only required for intra oral scanning.
2. Remove cover screw, healing abutment or similar from the implant or abutment upon which the Scan Body is intended to be placed.
3. Place the Scan Body onto the implant/ abutment in the patient's mouth or onto the analog when scanning a model.
4. Tighten the integrated screw slightly by hand using the Scan Body Driver.
5. Conduct the scan according to the scanner manufacturer's instructions to obtain a digital 3D model of the patient's mouth.
6. After scanning, the Scan Body is removed using the Scan Body Driver.

14.2 Design

1. Open the obtained digital 3D model in the scanners design software.
2. Select the relevant digital library which depend on the prosthetic solution used in the specific case.
3. Choose implant platform from the library corresponding to the scanned implant/analog or abutment.
4. Design the prosthetic restoration according to the design software's instructions. Be aware of that each specific prosthetic solution might have restricted design guidelines in accordance to maximum angulation and minimum wall thickness.

15 Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

16 Validity

Upon publication of this instruction for use, all previous versions are superseded.

17 Storage and Handling:

Recommended storage temperature range is 10-40°C (50-104°F) and the recommended relative humidity range is 0-80%

18 Disposal:

When disposed of the Scan Body may present an infection or microbial hazard and should therefore be disposed of as biological waste, according to national requirements.

19 Serious Incidents:

If any serious incidents occur in relation with this device, they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the

country where the incident occurred. A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat. Contact Elos Medtech Pinol A/S at: complaint.empi@elosmedtech.com

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

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complaint.empi@elosmedtech.com adresinden Elos Medtech Pinol A/S ile temasa geçin