Elos Accurate® Hybrid Base

Instruction for Use - Návod k použití –
Brugsanvisning – Gebruiksaanwijzing –
Kasutusjuhend – Käyttöohjeet – Mode
d'emploi – Gebrauchsanleitung – Οδηγίες
χρήσης – Használati útmutató – Istruzioni
per l'uso – Lietošanas norādījumi –
Naudojimo instrukcija – Bruksanvisning –
instrukcija Instrukcja użytkowania –
Instrucciones de uso

List of content

English – Instruction for Use	3
Čeština – Návod k použití	5
Dansk – Brugsanvisning	7
Nederlands – Gebruiksaanwijzing.....	9
Eesti keel – Kasutusjuhend.....	11
Suomi – Käyttöohjeet	13
Français — Mode d'emploi	15
Deutsch – Gebrauchsanleitung.....	17
Magyar – Használati útmutató	21
Italiano – Istruzioni per l'uso	23
Latviski — Lietošanas norādījumi	25
Lietuvių k. Naudojimo instrukcija.....	27
Norsk – Bruksanvisning	29
Język polski — Instrukcja użytkowania	31
Português — Instruções de uso (Portugal)	33
Español: Instrucciones de uso	35
Svenska – Bruksanvisning.....	37

English – Instruction for Use

Indications for Use: The Elos Accurate Hybrid Base Kit is intended for prosthetic restorations. The Elos Accurate Hybrid Base is used as an interface between a dental implant and a ceramic crown/bridge and will be attached to the implant using the included Prosthetic Screw and attached to the ceramic by cementing.

Device Description: The Elos Accurate Hybrid Base and Prosthetic Screw are manufactured from biocompatible titanium alloy grade 5 ELI (TiAl6V4 ELI). For specific product descriptions, please refer to individual product labels.

Contraindications: The Elos Accurate Hybrid Base is not intended for restorations requiring an angle more than 20° to the axis of the implant.

Allergies to the alloy, or contents of the alloy may very rarely occur.

Warnings: The connection to the implant system must never be changed or modified. The Elos Accurate Hybrid Base and included Prosthetic Screw are for single use only. Reuse can result in loss of functionality and/ or infections. Since the Elos Accurate Hybrid Base and Prosthetic Screw are small they must be handled with caution to avoid the patient swallowing or inhaling them.

When mounting the Elos Prosthetic Screw, it is important to use a manual screwdriver before using any kind of torque wrench.

Precautions: The Elos Accurate Hybrid Base must only be used and handled by dental professionals who have approved training in dentistry. The use of torque value higher than the recommended for the implant system according to the manufacturer's instruction may result in damage to the Elos Accurate Hybrid Base, the Prosthetic Screw and/or the implant. The use of torque values lower than those recommended for the implant system according to the manufacturer's instruction may result in loosening of the Elos Accurate Hybrid Base and/or the Prosthetic Screw. For any lab work and any try-in of the Elos Accurate Hybrid Base, a lab screw must be used. The included Prosthetic Screw may only be used to permanently fix the finished restoration.

MRI Statement: The Elos Accurate Hybrid Base and Prosthetic Screw have not been evaluated for safety, heating, migration or compatibility in the Magnetic Resonance Imaging (MRI) environment.

Potential Adverse Events: Potential adverse events associated with the use of Elos Accurate Hybrid Base products may include: failure to integrate, loss of integration and infection.

Note for Dental Professionals: Recommended torque value stated on label. See symbol .

Constructing the Custom Ceramic Abutment: The custom ceramic abutment must be constructed using approved design software and Elos Accurate Scanning Abutments. Internal geometry of the custom abutment is defined by the library file issued by Elos Medtech. For a ceramic abutment, minimum wall thickness is to be met according to the ceramic supplier's recommendation.

Milling a Custom Ceramic Abutment: The custom ceramic abutment must be manufactured with milling equipment adequate for the designed geometry of the ceramic abutment. A precise milling of the geometry requires millers of 1 mm in diameter or smaller.

Sintering the Custom Ceramic Abutment: The custom ceramic abutment should be sintered in accordance with the instruction issued by the ceramic material supplier.

Preparing the Custom Ceramic Abutment for Cementing: If the custom ceramic abutment is to be veneered, this must be done in accordance with the instruction issued by the ceramic material supplier. Internal geometry must be blasted in accordance with the instruction issued by ceramic material supplier. Before cementing, the surface should be thoroughly cleaned with alcohol or steam and dried.

Preparing the Elos Accurate Hybrid Base for Cementing: Blasting of Elos Accurate Hybrid Base is not necessary but optional. Before cementing, the surface should be thoroughly cleaned with alcohol or steam and dried. If the Elos Accurate Hybrid Base is to be blasted, it should be done in accordance with the instructions issued by the ceramic material supplier. In order to protect the connection to the implant, hold the Elos Accurate Hybrid Base in a lab / implant analog during blasting.

Cementing: Follow cementing procedure in accordance with instruction issued by cement supplier. In order to cement the abutment to the Elos Accurate Hybrid Base, insert the base into the lab / implant analog and fix it in place with the lab screw. Apply the cement in a circular fashion to the top edge of the contact surface of the Elos Accurate Hybrid Base and to the internal geometry of the abutment. Rotate the abutment gently as you push it over the Elos Accurate Hybrid Base until resistance is felt. Then turn it until its final position is reached. The custom abutment must fit the Elos Accurate Hybrid Base without gaps and voids. Any traces of cement should be removed immediately.


Polishing: Once the cement has set, the excess is removed with rubber polishing tools. The excess cement inside the screw channel must be carefully removed, ideally with the aid of a microscope. When polishing, fix the ceramic abutment in a lab / implant analog.

Cleaning and Sterilization: Prior to attachment of the Elos Accurate Hybrid Base /ceramic abutment and Prosthetic Screw to the implant, the Elos Accurate Hybrid Base / Abutment assembly and Prosthetic Screw must first be cleaned with water and brush followed by disinfection and sterilization. The recommended sterilization procedure is full cycle pre-vacuum steam sterilization at a temperature of 132 °C (270°F) for 4 minutes, 4 pulse minimum or full cycle gravity steam sterilization at 132°C (270°F) for 15 minutes minimum. The products are ready for use 20 minutes after autoclaving process.

Caution: US federal law restricts this product for sale by or on order of a dentist or physician. Elos Accurate Hybrid Bases may only be used by dentists/ physicians who have had the appropriate education and training.

Storage and Handling: Elos Accurate Hybrid Base Kit should be stored at room temperature.

Disposition: Must be disposed as biological waste.

	Catalogue number		Do not use if package is damaged
	Batch code		Prescription only
	Manufacturer		Date of manufacture
	Consult instructions for use		Medical Device
	Do not re-use		Recommended torque
	Non-sterile		

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