Elos Accurate® Hybrid Base™

Instruction for use

List of content

English - Instruction for use (English)	3
Български език – Инструкции за употреба (Bulgarian)	5
Hrvatski – Upute za upotrebu (Croatian)	8
Čeština - Návod k použití (Czech)	10
Dansk – Brugsanvisning (Danish)	12
Nederlands – Gebruiksaanwijzing (Dutch)	14
Eesti keel – Kasutusjuhised (Estonian)	17
Suomi – Käyttöohjeet (Finnish)	19
Français - Mode d'emploi (French)	21
Deutsch - Gebrauchsanleitung (German)	23
Ελληνικά - Οδηγίες χρήσης (Greek)	25
Magyar - Használati útmutató (Hungarian)	28
Íslenska – Notkunarleiðbeiningar (Icelandic)	30
Italiano - Istruzioni per l'uso (Italian)	32
Lietuvių k. - Naudojimo instrukcija (Lathuanian)	34
Latviski - Lietošanas norādījumi (Latvian)	36
Norsk – Bruksanvisning (Norwegian)	38
Polski - Instrukcja użytkowania (Polish)	40
Português - Instruções de uso (Portuguese)	43
Română – Instrucțiuni de utilizare (Romanian)	45
Slovenský jazyk – Návod na použitie (Slovak)	47
Slovenščina – Navodila za uporabo (Slovenian)	49
Español - Instrucciones de uso (Spanish)	51
Svenska - Anvisningar för användningen (Swedish)	53
Türkçe — Kullanma Talimatı (Turkish)	55

English - Instruction for use (English)

Elos Accurate® Hybrid Base™

1. Intended use

The Elos Accurate® Hybrid Base™ is intended for dental prosthetic restorations. The Hybrid Base™ is used as an interface between a dental implant or dental abutment and a dental restoration and will be attached to the implant or abutment using a prosthetic screw and attached to the dental restoration by cementing.

2. Product Description

The product consists of the Hybrid Base™ which is manufactured from biocompatible titanium alloy grade 5 ELI (TiAl₆V₄ ELI). The product is available for a variety of implant platforms and sizes. For specific product descriptions, please refer to individual product labels.

3. Indications

The Hybrid Base™ is directly or indirectly connected to the endosseous dental implant and is intended for use as an aid in prosthetic rehabilitations. The dental restoration can be cemented on the Hybrid Base™. A temporary restoration can be used prior to the insertion of the final restoration to maintain, stabilize and form the soft tissue during the healing phase. The temporary restoration may not be placed into occlusion. The Hybrid Base™ may be placed into occlusion when the implant is fully osseointegrated. The Hybrid Base™ Engaging is intended for single tooth dental restorations and the Hybrid Base™ Non-Engaging is intended for multiple tooth dental restorations.

4. Contraindications

- The Elos Accurate® Hybrid Base™ Non-Engaging is not intended for implants having a divergence angle more than 30° relatively to each other.
- The Elos Accurate® Hybrid Base™ Non-Engaging with a straight chimney (HBN-XXXXXH-1) is not intended for intra-oral cementation.
- The Elos Accurate® Hybrid Base™ Engaging is not intended for ceramic restorations angled more than 20° relative to the implant.

5. Warnings and precautions

- The Elos Accurate® Hybrid Base™ components must never be changed or modified.
- The component included in the product is for single-use only.
- Reuse of the component can result in loss of functionality and/or infections.
- The Elos Accurate® Hybrid Base™ must be attached to the implant or abutment using an Elos Prosthetic Screw or Elos Accurate® Hexalobular Prosthetic Screw.
- Tightening and loosening of the Prosthetic screw must be done by use of an appropriate screwdriver. Elos Accurate® Hexalobular Prosthetic Screws must be used with an Elos Accurate® Prosthetic Screwdriver 18mm, 26mm or 34mm (Ref. No. PS-AH18-1, PS-AH26-1 or PS-AH34-1).
- Since the Hybrid Base™ and Prosthetic Screw are small they must be handled with caution to avoid the patient from swallowing or inhaling them.
- The Hybrid Base™ and Prosthetic Screw have not been evaluated for safety, heating, migration or compatibility in the Magnetic Resonance Imaging (MRI) environment.
- Place implant-borne restorations only in occlusion when the implant is completely osseointegrated.
- Always place temporary restorations out of occlusion.
- Allergies to the titanium alloy grade 5 ELI (TiAl₆V₄ ELI) or contents of the alloy may very rarely occur.
- The Elos Accurate® Hybrid Base™ components must only be used and handled by dental professionals.
- The use of torque value higher or lower than the recommended for the implant system according to the manufacturer's instruction may result in damage to the Hybrid Base™, the Prosthetic screw and/or the implant
- Elos Accurate® Hybrid Base™ Non-Engaging is not intended for single tooth dental restorations.

6. Cleaning and disinfection

The Elos Accurate® Hybrid Base™ prosthetic components are delivered non-sterile. Prior to attachment of the dental restoration in the patient's mouth, it must be cleaned and disinfected followed by sterilization. Following cleaning and disinfection method is recommended:

1. Clean rinsing under flowing water while brushing outer and inner side with adequate brushes.
2. The pre-treated restoration can be cleaned by use of compatible cleaning solution.
3. The cleaned restoration can be disinfected in compatible disinfectant.
4. If the cleaned and disinfected dental restoration has to be packed, it should be packed in adequate packaging suitable for steam sterilization.

Cleaning and disinfection should be processed as specified by the manufacture. Cleaned and disinfected parts should be packed in adequate packaging according to instruction for use supplied by the dental restoration material manufacturer.

7. Sterilization

The Elos Accurate® Hybrid Base™ components are delivered non-sterile. The recommended sterilization procedure is full cycle pre-vacuum steam sterilization at a temperature of 132 °C (270°F) for 4 minutes or full cycle gravity steam sterilization at 132°C (270°F) for 15 minutes minimum. The products are ready for use 20 minutes after autoclaving process.

Caution: The device must be used immediately after sterilization. Do not store sterilized devices.

8. Procedure

8.1. Use and handling by dental technician

Gently open the outer packaging and take out the blister box. Gently open the blister to access the Hybrid Base™.

8.1.1. Manufacture the dental restoration

Recommended dental restoration material for Elos Accurate® Hybrid Base™ Engaging (permanent restorations):

- Zirconia.

Recommended dental restoration material for Elos Accurate® Hybrid Base™ Non-Engaging (permanent restorations):

- Zirconia.
- Cobalt-Chrome.
- Titanium (TiAl₆V₄ ELI).

Compatible design software: 3Shape, Exocad, and Dental Wings.

Caution:

- Ensure to follow the selected dental restoration material minimum wall thickness and parameters according to the material supplier's instruction for use.
- The dental restoration must be manufactured with milling equipment adequate for the designed geometry.
- A precise milling of the geometry requires millers of 1 mm in diameter or smaller.
- The dental restoration should if applicable be sintered in accordance with the instruction issued by the ceramic material supplier.

Using digital workflow (intra oral scanning):

1. Import the digitized patient situation from the dental clinical intraoral scan to the design software by use of Elos Accurate® Scan Body in order to identify the position and orientation of the respective implant.
2. Import library file from Elos Medtech's webpage and select relevant implant platform from the library in order to facilitate precise design of inner geometry of the dental restoration interfacing the Hybrid Base™

3. Design dental restoration in the design software or scanned teeth set up validated from the clinic in the patient's mouth.
4. Export STL file to milling equipment. The STL file contains the required milling template for the inner geometry of the dental restoration.
5. Process and finalize the dental restoration in milling equipment according to the manufacturer's instruction for use.
6. Create a digital working model in the design software.
7. Export the STL file from the design software and send the STL file to your 3D printer or external 3D print provider.
8. Place an Elos Accurate® Model Analog for printed models in the 3D printed working model.

Using semi-digital workflow (desktop scanning):

1. Create working model and place an Elos Accurate® Model Analog appropriately in the model.
2. Place an Elos Accurate® Scan Body in the model analog.
3. By use of the Elos Accurate® Scan Body, scan the working model.
4. Remaining steps as for the digital workflow.

8.1.2. Preparing the dental restoration for cementing

1. Blast inner geometry of the milled dental restoration interfacing the Hybrid Base™ with aluminum oxide 50-150 µm and blasting pressure of 2 bar.
2. Before cementing, clean the surface thoroughly with alcohol.

Note: Blasting of Hybrid Base™ is not necessary, but if blasted make sure to protect the implant/abutment interface.

8.1.3. Cementing the dental restoration

1. Seal the screw channels with wax.
2. Apply self-adhesive dental cement on the the Hybrid Base™. Use only suitable self-adhesive cementation systems for the dental restoration material used. Follow the instruction for use of both the dental restoration material and cement material manufacturer. Bond the milled dental restoration to the Hybrid Base™ placed in the working model. The milled dental restoration must fit the Hybrid Base™ without gaps and voids.
3. Immediately remove excess cement from the Hybrid Base™.

9. Use and handling for the dentist

The dentist receives the final dental restoration / working model with the Prosthetic screw from the dental lab.

1. Identify and unpack the Prosthetic screw(s).
2. Clean, disinfect and sterilize dental restoration and Prosthetic screw(s) according to this instruction for use.
3. Remove the healing cap, closure screw or temporary restoration from patient's mouth.
4. Gently, insert the dental restoration into the patient's mouth in proper position to the implant(s) or abutment(s).
5. Place corresponding Prosthetic screw(s) in the dental restoration and tighten screws using torque values specified by the implant manufacturer.
6. Additionally, if cementing of the dental restoration and the Hybrid Base™ has to be done intraorally, place Hybrid Base™ in the milled dental restoration with correct rotation.
7. Seal the screw channels with wax.
8. Apply self-adhesive dental cement on the the Hybrid Base™ Use only suitable self-adhesive cementation systems for the dental restoration material used. Follow the instruction for use of both the dentalrestoration material and cement material manufacturer. Recommended cement is described in this instruction for use.
9. Bond the milled dental restoration to the Hybrid Base™. The milled dental restoration must fit the Hybrid Base™ without gaps and voids. Immediately remove excess cement from the Hybrid Base™.

10. Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

11. Validity

Upon publication of this instruction for use, all previous versions are superseded.


12. Storage and Handling

Elos Accurate® Hybrid Base™ should be stored at room temperature.

13. Disposal

The dental restoration must be disposed as biological waste.

14. Symbols

	Catalogue number
	Batch code
	Manufacturer
	Consult instructions for use
	Do not re-use
	Prescription only
	Do not use if package is damaged
	Non sterile
	Recommended torque
	Date of manufacture
	Medical Device

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