Universal Base Non-Engaging

Instruction for use

Powered by Elos Accurate®

List of content

English – Instruction for use (English)	3
Български – Инструкция за употреба (Bulgarian)	6
Hrvatski – upute za uporabu (Croatian)	9
Čeština – Návod k použití (Czech)	12
Dansk – Brugsanvisning (Danish)	15
Nederlands – Gebruiksaanwijzing (Dutch)	18
Eesti keel – Kasutusjuhend (Estonian)	21
Suomi – Käyttöohjeet (Finnish)	24
Français – Mode d'emploi (French)	27
Deutsch – Gebrauchsanleitung (German)	30
Ελληνικά – Οδηγίες χρήσης (Greek)	33
Magyar – Használati útmutató (Hungarian)	36
Íslenska – notkunarleiðbeiningar (Icelandic)	39
Italiano – Istruzioni per l'uso (Italian)	42
Lietuvių k. – naudojimo instrukcija (Lithuanian)	45
Latviski – Lietošanas norādījumi (Latvian)	48
Norsk – Bruksanvisning (Norwegian)	51
Polski – Instrukcja użytkowania (Polish)	54
Português – Instruções de uso (Portuguese)	57
Română – Instrucțiuni de utilizare (Romanian)	60
Slovenčina – Návod na použitie (Slovak)	63
Slovenščina – navodila za uporabo (Slovenian)	66
Español – Instrucciones de uso (Spanish)	69
Svenska – Bruksanvisning (Swedish)	72
Türkçe – Kullanma Talimatı (Turkish)	75

English – Instruction for use (English) Universal Base Non-Engaging

1. Intended purpose

The Universal Base is intended for dental prosthetic restorations. The Universal Base is used as an interface between an endosseous dental implant and a dental restoration and will be attached to the implant using a prosthetic screw and attached to the dental restoration by cementing.

Universal Base Non-engaging are applicable for multiple tooth restorations.

2. Product description

The product consists of the Universal Base and prosthetic screw which is manufactured from biocompatible titanium alloy (TiAl₆V₄ ELI). Depending on the specific variant, the Prosthetic Screw includes biocompatible coatings.

The product is available for a variety of implant platforms and sizes. For specific product descriptions, please refer to individual product labels.

3. Chemical composition

Chemical composition of materials and coatings used for the Universal Base and prosthetic screw is presented in table 1-3

Table 1: Chemical composition of Universal Base prosthetic screws made of anodized titanium alloy

Material	Composition	Patient contacting	CAS nr.
Anodized Titanium alloy (TIODIZE Type II)	33.3 at.% Ti + 66.6 at.% O (TiO ₂)	Yes (Outer layer)	1317-80-2
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	No (Inner layer)	99906-66-8

Table 2: Chemical composition of Universal Base prosthetic screws made of titanium alloy with DLC coating

Material	Composition	Patient contacting	CAS nr.
Tungsten carbide and Amorphous carbon (Diamond like carbon)	80 at.% C + 16.8 at.% W with residual Ni	Yes (Outer layer)	12070-12-1 (WC) + 7782-42-5 (C)
Tungsten carbide binding layer	50 at.% W + 50 at.% C (WC)	No (Inner layer)	12070-12-1
Chrome (II) binding layer	100 wt.% Cr (II)	No (Inner layer)	7440-47-3
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	No (Inner layer)	99906-66-8

Table 3: Chemical composition of the Universal Base Abutment & prosthetic screws made of titanium alloy without a coating

Material	Composition	Patient contacting	CAS nr.
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	Yes (Pure metal)	99906-66-8

4. Basic UDI:

Universal Base Non-engaging: 5712821011001UD

5. Indication for use

The Universal Base is indicated for patients with edentulism, who need single or multiple dental restorations attached to a dental implant.

6. Intended patient population

The Universal Base is intended to be used in patients that are considered eligible for dental implant treatment. No special training of the patient is required.

7. Intended users

For the preparation of the dental restoration the Universal Base should only be handled by dental technicians who are experienced within prosthetics for dental implantology.

For the installation of the dental restoration in patients the Universal Base should only be used by dentists who are experienced within dental implantology.

7.1. Use environment

Dental labs need to have 3D milling equipment suitable for milling of customized dental restorations.

8. Clinical benefits

The Universal Base can successfully be used to establish a prosthetic restoration to treat patients with edentulism.

9. Expected lifetime

Expected lifetime of up to 25 years, corresponding to the state of the art as assessed in the clinical evaluation.

10. Contraindications

- The Universal Base Non-Engaging is not intended for implants having a divergence angle more than 30° relatively to each other.

11. Warnings and precautions

- The Universal Base components must never be changed or modified.
- The component included in the product is for single-use only.
- Reuse of the component can result in loss of functionality and/or infections.
- The Universal Base must be attached to the implant or abutment using a Prosthetic Screw.
- Tightening and loosening of the Prosthetic screw must be done by use of an appropriate screwdriver.
- Since the Universal Base and Prosthetic Screw are small they must be handled with caution to avoid the patient from swallowing or inhaling them.
- Place implant-borne restorations only in occlusion when the implant is completely osseointegrated.
- Always place temporary restorations out of occlusion.
- Allergies to the titanium alloy grade 5 ELI (TiAl₆V₄ ELI) or contents of the alloy may very rarely occur.
- The Universal Base components must only be used and handled by dental professionals.
- The use of torque value higher or lower than the recommended for the implant system according to the manufacturer's instruction may result in damage to the Universal Base, the Prosthetic screw and/or the implant
- Universal Base Non-Engaging is not intended for single tooth dental restorations.
- When mounting the Prosthetic Screw, it is important to use a manual screwdriver before using any kind of torque wrench.
- When intra oral cementation is intended, make sure that insertion direction of the restoration allows it.

12. Compatibility information

The Universal Base is available for a variety of implant platforms and sizes. The individual product label refers to which implant platform and size the product is compatible with.

13. Cleaning and Sterilization

The Universal Base prosthetic components are delivered non-sterile. Prior to attachment of the dental restoration in the patient's mouth, it must be cleaned followed by sterilization. Further instruction can be found in "Cleaning and sterilization guideline" at:

<https://elosmedtech.com/IFU/>

- Clean according to the guidance in the "Cleaning and sterilization guideline".
- Sterilize according to the guidance in the "Cleaning and sterilization guideline". Summarized below.

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Caution: The device must be used immediately after sterilization. Do not store sterilized devices.

14. Procedure

14.1. Use and handling by dental technician

Gently open the outer packaging and take out the blister box. Gently open the blister to access the Universal Base.

14.1.1. Manufacture the dental restoration

Recommended dental restoration material for Universal Base Non-Engaging (permanent restorations):

- Zirconia.
- Titanium (TiAl₆V₄ ELI).

Compatible design software: 3Shape, Dental Wings, Exocad and DTX Studio Lab.

Caution:

- Ensure to follow the selected dental restoration material minimum wall thickness and parameters according to the material supplier's instruction for use.
- The dental restoration must be manufactured with milling equipment adequate for the designed geometry.
- A precise milling of the geometry requires millers of 1 mm in diameter or smaller.
- The dental restoration should if applicable be sintered in accordance with the instruction issued by the ceramic material supplier.

Using digital workflow (intra oral scanning):

1. Import the digitized patient situation from the dental clinical intraoral scan to the design software by use of Elos Accurate® Scan Body in order to identify the position and orientation of the respective implant.
2. For DTX Studio Lab: the implant libraries are automatically included in the software installer.
3. For 3Shape, Dental Wings and Exocad: download libraries from www.nobelbiocare.com/international/en/local-production.
4. Select relevant implant platform from the library in order to facilitate precise design of inner geometry of the dental restoration interfacing the Universal Base.
5. Design dental restoration in the design software or scanned teeth set up validated from the clinic in the patient's mouth.
6. Export STL file to milling equipment. The STL file contains the required milling template for the inner geometry of the dental restoration.
7. Process and finalize the dental restoration in milling equipment according to the manufacturer's instruction for use.
8. Create a digital working model in the design software.
9. Export the STL file from the design software and send the STL file to your 3D printer or external 3D print provider.
10. Place an Elos Accurate® Model Analog for printed models in the 3D printed working model.

Using semi-digital workflow (desktop scanning):

1. Create working model and place an Elos Accurate® Model Analog appropriately in the model.
2. Place an Elos Accurate® Scan Body in the model analog.
3. By use of the Elos Accurate® Scan Body, scan the working model.
4. Remaining steps as for the digital workflow.

14.1.2. Preparing the dental restoration for cementing

1. Blast inner geometry of the milled dental restoration interfacing the Universal Base with aluminum oxide 50-150 µm and blasting pressure of 2 bar.
2. Before cementing, clean the surface thoroughly with alcohol.

Note: Do not blast the Universal Base.

14.1.3. Cementing the dental restoration

1. Seal the screw channels with wax.
2. Apply a primer to the bonding surface of zirconia and Universal Base and let it react according to manufactures instruction.
3. Apply self-adhesive dental cement on the the Universal Base. Follow the instruction for use of both the dental restoration material and cement material manufacturer. Bond the milled dental restoration to the Universal Base placed in the working model or in the patients mouth. The milled dental restoration must fit the Universal Base without gaps and voids.
4. Immediately remove excess cement from the Universal Base.

One of following bonding agents is recommended.

Primer / Cement	Restoration material	Intra oral cementation
Kuraray Noritake CLEARFIL™ CERAMIC PRIMER PLUS / PANA VIA™ V5	Zirconia	Yes
Ivoclar Vivadent® Monobond Plus / Multilink® Hybrid Abutment	Zirconia Titanium	No
3M™ Scotchbond™ Universal Plus Adhesive / RelyX™ Universal Resin Cement	Zirconia	Yes

Further instruction can be found in "Cementing guideline" at <https://elosdental.com/training-and-support/guides-our-products>

15. Use and handling for the dentist

The dentist receives the final dental restoration / working model with the Prosthetic screw from the dental lab.

1. Identify and unpack the Prosthetic screw(s).
2. Clean, disinfect and sterilize dental restoration and Prosthetic screw(s) according to this instruction for use.
3. Remove the healing cap, closure screw or temporary restoration from patient's mouth.
4. Gently, insert the dental restoration into the patient's mouth in proper position to the implant(s) or abutment(s).
5. Place corresponding Prosthetic screw(s) in the dental restoration and tighten screws using torque values specified by the implant manufacturer.
6. Seal the screw channels with wax.
7. Apply self-adhesive dental cement on the the Universal Base Use only suitable self-adhesive cementation systems for the dental restoration material used. Follow the instruction for use of both the dental restoration material and cement material manufacturer. Recommended cement is described in this instruction for use.
8. Bond the milled dental restoration to the Universal Base. The milled dental restoration must fit the Universal Base without gaps and voids. Immediately remove excess cement from the Universal Base.

16. Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

17. Validity

Upon publication of this instruction for use, all previous versions are superseded.

18. Storage and Handling

Universal Base should be stored at room temperature.

19. MRI Safety Information

<p>MR Safety Information</p> <p>A person with an Universal Base Non-Engaging Abutment, associated dental implant and prosthetic screw may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.</p>	
Device Name	Universal Base Non-Engaging
Static Magnetic Field Strength (B ₀)	1.5 T or 3.0 T
Maximum Spatial Field Gradient	20 T/m (2,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil
Operation Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	2 W/kg whole-body average SAR for 15 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR Image Artifact	The presence of passive implant devices may produce an image artifact that scales with the device size

	Non sterile
	Recommended torque
	Date of manufacture
	Medical Device
	Unique Device Identifier
	Keep dry
	Keep away from sunlight
	MR Conditional

20. Disposal

The dental restoration must be disposed as biological waste.

21. Serious Incidents:

If any serious incidents occur in relation with this device, they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the country where the incident occurred.

A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat.

Contact Elos Medtech Pinol A/S at:

complaint.empi@elosmedtech.com

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. URL: <https://ec.europa.eu/tools/eudamed>.

22. Symbols

	Catalogue number
	Batch code
	Manufacturer
	Consult instructions for use
	Do not re-use
	Prescription only
	Do not use if package is damaged

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Şu adres üzerinden Elos Medtech Pinol A/S ile iletişime geçin: