Cleaning- and sterilization guideline for Elos Medtech reusable medical devices

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Disclaimer of liability: Elos Medtech Pinol products may only be used according to the instructions and recommendations. The user of Elos Medtech Pinol's products should determine their suitability for particular patients and indications. Elos Medtech Pinol disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or damages arising out of or in connection with any errors in professional judgement or practice in the use or placement of the Elos Medtech Pinol products.

Please note that not all products are regulatory cleared/released for sale in all markets even though the language is present in this instruction.

Elos Medtech Pinol A/S

Engvej 33

DK - 3330 Gørlev

www.elosmedtech.com

This electronic instruction for use can be accessed at <https://elosmedtech.com/IFU/>

Document ID: DEV-01653 version 3.0 released 2020/03/16

Cleaning- and sterilization guideline for Elos Medtech reusable medical devices – (English)

Content

1	Purpose	3
2	Scope	3
3	Initial treatment at the point of use	3
4	Preparation for cleaning	3
5	Cleaning	3
6	Automated disinfection.....	4
7	Drying	4
8	Maintenance, Inspection and Testing.....	4
9	Packaging.....	4
10	Sterilization.....	4
11	Storage	5
12	Manufacturer contact.....	5

1 Purpose

The purpose of guideline is to describe Elos Medtech's cleaning, sterilization and dry time procedure for Elos Medtech's reusable medical devices.

2 Scope

This guideline specifies the recommended methods to be used for safe and effective manual- and automated cleaning and sterilization. The guideline was made in compliance to international regulatory standards and guidance documents as listed below.

- AAMI TIR12 (2010)
- AAMI TIR30 (2011/(R) 2016)
- FDA Guidance (March 17, 2015, updated June 09, 2017)
- ISO 17664 (2017)
- ANSI/AAMI ST79 (2017)

3 Initial treatment at the point of use

- To prevent any contamination from drying it is recommended that cleaning- and sterilization takes place not later than 2 hours after use.
- Keep instruments moist after use to prevent soil from drying on them.
- Contaminated instruments should be transported to the area for cleaning in a clean and closed box/container that avoids spreading of the contamination.

4 Preparation for cleaning

Where applicable, disassemble instruments until each device consists of one part.

5 Cleaning

- Prepare a solution of a mild alkaline detergent Neodisher MediClean Forte 0.5%, temperature 15-20 °C.
- Completely submerge the instruments in the cleaning solution for five (5) minutes. Whilst submerged, brush the instrument with attention to the critical and hard-to reach areas with a soft-bristled brush for a minimum of 30 seconds and flush hollow instruments using a syringe with 0.5 mL of cleaning solution.
- Rinse one (1) minute under cold running tap water (15-20°C)
- Place the devices into the washing racks without overloading them. Then place the racks on the wash supports of the washer-disinfector. Avoid any contact between the devices, which may cause damage during washing
- Using a valid washer disinfector (according to ISO 15883 requirements) and an alkaline enzymatic detergent such as Neodisher MediClean Forte, use the minimal cycle parameters set points mentioned in the following table.

Table 1 Procedure for cleaning

Cycle	Minimum Time	Minimum Temperature	Type of Detergent/Water
Pre-Cleaning	2 minutes	Cold (< 45°C)	Tap water
Wash	10 minutes	Heated (50-60°C)	Neodisher MediClean Forte (0.2-1%)
Rinse	2 minutes	Cold (< 45°C)	Critical water ¹
Thermal Rinse	5 minutes	Heated (90°C)	Critical water ¹
Dry	25 minutes	Heated (90°C)	Not applicable

- At the end of the program remove the instruments.
- If not completely dry, dry instruments with single use, non-linting wipe

6 Automated disinfection

For Europe, the thermal disinfection corresponds to the Thermal Rinse phase of five (5) minutes at a minimum of 90°C, listed in Table 2

7 Drying

If not completely dry, dry instruments with single use, non-linting wipe or a towel or by air gun with compressed air.

8 Maintenance, Inspection and Testing

Any instrument showing visible damage or wear must be exchanged. End of life of instrument is normally determined by wear and damage due to use. End of life of instruments is exceeded if visual inspection shows unacceptable deterioration like wear, corrosion, discoloration, pitting crackles etc.

9 Packaging

- Instruments should be sterilized disassembled according to their Instructions For Use
- Place the instruments in the double sterilization packaging (pouches made of paper or plastic film, conforming to ISO 13060, DIN EN 258 and ANSI/AAMI ST79
- For USA use FDA-cleared pouches
- Physical-chemical indicators may be used on or in the packaging system
- For USA use FDA-approved wraps

10 Sterilization

- Use a steam sterilizer conforming to EN 13060, DIN 285, and /or ST79
- Use the sterilization parameters provided in the table below, which have been validated

¹ As per AAMI TIR34, water extensively treated usually by a multistep treatment process that may include a carbon bed, softening, DI, and RO or distillation, to ensure that the microorganisms and the inorganic and organic material are removed from the water

Note: For US: Steam Sterilization Cycle 132°C (270°F)/4min.

Table 2 Dynamic-air-removal Steam Sterilization Cycle, Exposure Time, Temperature & Drying Time

Procedure	Dynamic-air-removal-sterilization Cycle	
Exposure time	4 min.	3 min.
Temperature	132°C (270°F)	134°C (273°F)
Drying time	20 min.	

11 Storage

- Keep devices in sterilization packaging in a dry and clean environment.
- Sterility cannot be guaranteed if packaging becomes open, damaged or wet.
- Check the packaging and the medical device before using them (package integrity, no excessive humidity and validity period)

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