# CORiTEC® Prosthetic Screw

Instruction for use - Инструкция за употреба  
- Návod k použití - Brugsanvisning -  
Gebrauchsanleitung - Kasutusjuhend -  
Οδηγίες χρήσης - Instrucciones de uso -  
Mode d'emploi - Upute za upotrebu -  
Notkunarleiðbeiningar - Istruzioni per l'uso -  
Lietošanas norādījumi - Naudojimo instrukcija  
- Használati útmutató - Gebruiksaanwijzing -  
Bruksanvisning - Instrukcja użytkowania -  
Instruções de uso - Instrucțiuni de utilizare -  
Návod na použitie - Navodila za uporabo -  
Käyttöohjeet - Bruksanvisning - Kullanma  
Talimatı

## List of content

English – Instruction for use (English).....	3
Български – Инструкция за употреба (Bulgarian) .....	5
Hrvatski – upute za uporabu (Croatian) .....	8
Čeština – Návod k použití (Czech).....	10
Dansk – Brugsanvisning (Danish).....	12
Nederlands – Gebruiksaanwijzing (Dutch) .....	14
Eesti keel – Kasutusjuhend (Estonian).....	16
Suomi – Käyttöohjeet (Finnish).....	18
Français – Mode d’emploi (French) .....	20
Deutsch – Gebrauchsanleitung (German).....	23
Ελληνικά – Οδηγίες χρήσης (Greek) .....	26
Magyar – Használati útmutató (Hungarian).....	29
Íslenska – notkunarleiðbeiningar (Icelandic) .....	31
Italiano – Istruzioni per l’uso (Italian).....	33
Lietuvių k. – naudojimo instrukcija (Lithuanian).....	36
Latviski – Lietošanas norādījumi (Latvian) .....	38
Norsk – Bruksanvisning (Norwegian).....	40
Polski – Instrukcja użytkowania (Polish) .....	42
Português – Instruções de uso (Portuguese).....	45
Română – Instrucțiuni de utilizare (Romanian) .....	47
Slovenčina – Návod na použitie (Slovak).....	49
Slovenščina – navodila za uporabo (Slovenian).....	51
Español – Instrucciones de uso (Spanish) .....	53
Svenska – Bruksanvisning (Swedish) .....	55
Türkçe – Kullanma Talimatı (Turkish) .....	57

**English – Instruction for use (English)**  
**CORiTEC® Prosthetic Screw**

**1. Intended purpose**

The CORiTEC® Prosthetic Screw is intended for fixation of prosthetic restorations to endosseous dental implants.

**2. Product description**

The product consists of the Prosthetic Screw which is manufactured from biocompatible titanium (Grade 4) or titanium alloy (TiAl<sub>6</sub>V<sub>4</sub> ELI). Depending on the specific variant, the CORiTEC® Prosthetic Screw includes biocompatible coatings.

The product is available for a variety of implant platforms and sizes. For specific product descriptions, please refer to individual product labels.

**3. Chemical composition**

Chemical composition of materials and coatings used for the CORiTEC® Prosthetic Screw is presented in table 1-4.

*Table 1: Chemical composition of prosthetic screws made of titanium alloy with MediCarb coating*

Material	Composition	Patient contacting	CAS nr.
Amorphous carbon (Diamond like carbon)	100 at. % C	Yes (Outer layer)	7782-42-5
Chromium nitride (CrN)	50 at. % Cr (III) and 50 at. % N	No (Inner layer)	24094-93-7
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	No (Inner layer)	99906-66-8

*Table 2: Chemical composition of prosthetic screws made of titanium alloy with MediGold coating*

Material	Composition	Patient contacting	CAS nr.
Gold	97.89 at. % Au with residual Ni	Yes (Outer layer)	7440-57-5
Titanium Nitride	50 at. % Ti and 50 at. % N	No (Inner layer)	25583-20-4
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	No (Inner layer)	99906-66-8

*Table 3: Chemical composition of prosthetic screws made of titanium alloy without a coating*

Material	Composition	Patient contacting	CAS nr.
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	Yes (Pure metal)	99906-66-8

*Table 4: Chemical composition of prosthetic screws made of titanium grade 4 without a coating*

Material	Composition	Patient contacting	CAS nr.
Titanium Grade 4	≥ 98.6 wt.% Ti, ~ 0.50 wt.% Fe and 0.4 wt.% O	Yes (Pure metal)	7440-32-6

**4. Basic UDI:**

Basic UDI: 5712821006001V5

**5. Indication for use**

The CORiTEC® Prosthetic Screw is indicated for patients with edentulism, who need single or multiple dental restorations attached to a dental implant.

**6. Intended patient population**

The CORiTEC® Prosthetic Screw is intended to be used on patients that are considered eligible for dental implant treatment. No special training of the patient is required.

**7. Intended users**

For the installation of the dental restoration in patients the CORiTEC® Prosthetic Screw should only be used by dentists who are experienced within dental implantology.

**8. Use environment**

Dental labs need to have 3D milling equipment suitable for milling of customized dental restorations.

**9. Clinical benefits**

The CORiTEC® Prosthetic Screws can successfully be used for fixation of prosthetic restorations to treat patients with edentulism.

**10. Expected life lifetime**

Expected lifetime of up to 25 years, corresponding to the state of the art as assessed in the clinical evaluation.

**11. Contraindications**

- The CORiTEC® Prosthetic Screw is not intended for any lab work or try-in of the prosthetic restoration.

**12. Warnings and precautions**

- The component included in the product is for single-use only.
- Reuse of the component can result in loss of functionality and/ or infections.
- Tightening and loosening of the Prosthetic screw must be done by use of an appropriate screwdriver.
- When mounting the Prosthetic Screw, it is important to use a manual screwdriver before using any kind of torque wrench.
- Since the Prosthetic Screw are small they must be handled with caution to avoid the patient from swallowing or inhaling them.
- Allergies to the titanium/ titanium alloy (TiAl<sub>6</sub>V<sub>4</sub> ELI) or contents of the alloy may very rarely occur.
- The CORiTEC® Prosthetic Screw components must only be used and handled by dental professionals.
- For any lab work and any try-in of the prosthetic restoration, a lab screw must be used. The CORiTEC® Prosthetic Screw may only be used to permanently fix the finished restoration
- The use of torque value higher or lower than the recommended for the implant system according to the manufacturer's instruction may result in damage to the prosthetic restoration, the prosthetic screw and/or the implant

**13. Compatibility information**

The Prosthetic Screw is available for a variety of implant platforms and sizes. The individual product label refers to which implant platform and size the product is compatible with.

**14. Cleaning and Sterilization**

The CORiTEC® Prosthetic Screw is delivered non-sterile. Prior to installation of the Prosthetic Screw in the patient's mouth, it must be cleaned followed by sterilization. Further instruction can be found in "Cleaning and sterilization guideline" at:

<https://elosmedtech.com/IFU/>

- Clean according to the guidance in the "Cleaning and sterilization guideline".
- Sterilize according to the guidance in the "Cleaning and sterilization guideline". Summarized below.

Procedure	Dynamic-air-removal sterilization Cycle
Exposure time	3 min.
Temperature	134° C (273°F)
Drying time	20 min.

**Caution:** The device must be used immediately after sterilization. Do not store sterilized devices.

Elos Medtech Pinol A/S  
Engvej 33  
DK - 3330 Gørlev  
[www.elosmedtech.com](http://www.elosmedtech.com)

**15. Procedure**

**Use and handling for the dentist**

Gently open the packaging to access the Prosthetic screw.

1. Gently open the packaging to access the Prosthetic screw(s).
2. Clean, disinfect and sterilize the Prosthetic screw(s) according to this instruction for use.
3. Place corresponding Prosthetic screw(s) in the dental restoration and tighten screws using torque values specified by the implant manufacturer.

**16. Further information**

For additional information about the use of CORiTEC® Medtech products, please contact your local sales representative.

**17. Validity**

Upon publication of this instruction for use, all previous versions are superseded.

**18. Storage and Handling**

CORiTEC® Prosthetic Screw should be stored at room temperature.

**19. MRI Safety Information**

 MRI Safety Information	
A person with an CORiTEC® prosthetic screw, associated dental implant and abutment may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
Device Name	CORiTEC® Prosthetic Screw
Static Magnetic Field Strength (B <sub>0</sub> )	1.5 T or 3.0 T
Maximum Spatial Field Gradient	20 T/m (2,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil
Operation Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR Image Artifact	The presence of passive implant devices may produce an image artifact that scales with the device size

**20. Disposal**

The CORiTEC® Prosthetic Screw must be disposed as biological waste.

**21. Serious Incidents:**

If any serious incidents occur in relation with this device they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the country where the incident occurred.

A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat.


Contact Elos Medtech Pinol A/S at:

[complaint.empi@elosmedtech.com](mailto:complaint.empi@elosmedtech.com)

The SSCP is available in the European database on medical

devices (Eudamed), where it is linked to the Basic UDI-DI. URL: <https://ec.europa.eu/tools/eudamed>.

**22. Symbols**

	Catalogue number
	Batch code
	Manufacturer
	Consult instructions for use
	Do not re-use
	Prescription only
	Do not use if package is damaged
	Non sterile
	Recommended torque
	Date of manufacture
	Medical Device
	Unique Device Identifier
	Keep dry
	Keep away from sunlight
	MR Conditional

Elos Medtech Pinol A/S  
Engvej 33  
DK - 3330 Gørølse  
[www.elosmedtech.com](http://www.elosmedtech.com)

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Şu adres üzerinden Elos Medtech Pinol A/S ile iletişime geçin: