CORiTEC® preMill

Instruction for use - Инструкция за
употреба – Návod k použití –
Brugsanvisning – Gebrauchsanleitung –
Kasutusjuhend – Instrucciones de uso –
Οδηγίες χρήσης – Mode d'emploi – Upute za
upotrebu – Istruzioni per l'uso – Lietošanas
norādījumi – Naudojimo instrukcija –
Használati útmutató – Notkunarleiðbeiningar
– Gebruiksaanwijzing – Bruksanvisning –
Instrukcja użytkowania – Instruções de uso –
Instrucțiuni de utilizare – Návod na použitie –
Navodila za uporabo – Käyttöohjeet –
Bruksanvisning – Kullanma Talimatı

List of content

English – Instruction for use (English)	3
Български – Инструкция за употреба (Bulgarian)	6
Hrvatski – upute za uporabu (Croatian).....	9
Čeština – Návod k použití (Czech).....	12
Dansk – Brugsanvisning (Danish)	15
Nederlands – Gebruiksaanwijzing (Dutch)	18
Eesti keel – Kasutusjuhend (Estonian).....	21
Suomi – Käyttöohjeet (Finnish)	24
Français – Mode d'emploi (French).....	27
Deutsch – Gebrauchsanleitung (German).....	30
Ελληνικά – Οδηγίες χρήσης (Greek)	33
Magyar – Használati útmutató (Hungarian).....	36
Íslenska – notkunarleiðbeiningar (Icelandic)	39
Italiano – Istruzioni per l'uso (Italian)	42
Lietuvių k. – naudojimo instrukcija (Lithuanian).....	45
Latviski – Lietošanas norādījumi (Latvian)	48
Norsk – Bruksanvisning (Norwegian)	51
Polski – Instrukcja użytkowania (Polish).....	54
Português – Instruções de uso (Portuguese)	57
Română – Instrucțiuni de utilizare (Romanian)	60
Slovenčina – Návod na použitie (Slovak).....	63
Slovenščina – navodila za uporabo (Slovenian).....	66
Español – Instrucciones de uso (Spanish)	69
Svenska – Bruksanvisning (Swedish)	72
Türkçe – Kullanma Talimatı (Turkish)	75

English – Instruction for use (English)

CORiTEC® preMill

1 Intended purpose

The CORiTEC® preMill is intended for dental prosthetic restorations. The CORiTEC® preMill is used as an interface between an endosseous dental implant and a dental restoration and will be attached to the implant using a prosthetic screw and attached to the dental restoration by cementing.

2 Product Description

The CORiTEC® preMill is intended for production of patient-specific abutments. The product consists of the CORiTEC® preMill which is manufactured from biocompatible titanium alloy (TiAl₆V₄ ELI).

The product is available for a variety of implant platforms and sizes. For specific product descriptions, please refer to individual product labels.

3 Chemical composition

Chemical composition of materials used for the CORiTEC® preMill is presented in table 1.

Table 1: Chemical composition of CORiTEC® preMill made of titanium alloy

Material	Composition	Patient contacting	CAS nr.
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	Yes (Pure metal)	99906-66-8

4 Basic UDI:

CORiTEC® preMill: 5712821005001UW

5 Indication for use

The CORiTEC® preMill is indicated for patients with edentulism, who need single dental restorations attached to a dental implant.

6 Intended patient population

The CORiTEC® preMill is intended to be used in patients that are considered eligible for dental implant treatment. No special training of the patient is required.

7 Intended users

For the preparation of the dental restoration the CORiTEC® preMill should only be handled by dental technicians who are experienced within prosthetics for dental implantology and milling of customized abutments from CORiTEC® preMills.

For the installation of the dental restoration in patients the CORiTEC® preMill should only be used by dentists who are experienced within dental implantology.

7.1 Use environment

Dental labs need to have 3D milling equipment suitable for milling of customized dental restorations.

8 Clinical benefits

The CORiTEC® preMill can successfully be used to establish a prosthetic restoration to treat patients with edentulism.

9 Expected lifetime

Expected lifetime of up to 25 years, corresponding to the state of the art as assessed in the clinical evaluation.

10 Contraindications

- The CORiTEC® preMill is not intended for restorations angled more than 30° relative to the implant.

11 Warnings and precautions

- The connection interface of the CORiTEC® preMill must never be changed or modified.
- The component included in the product is for single-use only.
- Reuse of the component can result in loss of functionality and/ or infections.
- The CORiTEC® preMill must be attached to the implant or abutment using an CORiTEC® Prosthetic Screw.
- Tightening and loosening of the Prosthetic screw must be done by use of an appropriate screwdriver.
- When mounting the Prosthetic Screw, it is important to use a manual screwdriver before using any kind of torque wrench.
- Since the CORiTEC® preMill and Prosthetic Screw are small they must be handled with caution to avoid the patient from swallowing or inhaling them.
- Place implant-borne restorations only in occlusion when the implant is completely osseointegrated.
- Always place temporary restorations out of occlusion.
- Allergies to the titanium alloy (TiAl₆V₄ ELI) or contents of the alloy may very rarely occur.
- The CORiTEC® preMill components must only be used and handled by dental professionals.
- The use of torque value higher or lower than the recommended for the implant system according to the manufacturer's instruction may result in damage to the CORiTEC® preMill, the Prosthetic screw and/or the implant
- Risk of fire when machining titanium without cutting emulsion
- Titanium chips should be disposed of in a suitable waste container (fire-proof)

12 Compatibility information

The CORiTEC® preMill is available for a variety of implant platforms and sizes. The individual product label refers to which implant platform and size the product is compatible with.

13 Cleaning and Sterilization

The final CORiTEC® components customized prosthetic components are delivered non-sterile. Prior to attachment of the dental restoration in the patient's mouth, it must be cleaned followed by sterilization. Further instruction can be found in "Cleaning and sterilization guideline" at:

<https://elosmedtech.com/IFU/>

- Clean according to the guidance in the "Cleaning and sterilization guideline".
- Sterilize according to the guidance in the "Cleaning and sterilization guideline". Summarized below.

Procedure	Dynamic-air-removal sterilization Cycle
Exposure time	3 min.
Temperature	134° C (273°F)
Drying time	20 min.

Caution: The device must be used immediately after sterilization. Do not store sterilized devices.

14 Procedure

14.1 Use and handling by dental technician

Gently open the packaging to access the CORiTEC® preMill.

14.1.1 Manufacture the customized abutment

Compatible design software: 3Shape, Exocad, and Dental Wings.

Caution:

- Ensure that wall thickness of customized abutment is minimum 0.4 mm.
- The screw head must not be exposed when designing the customized abutment.
- The customized abutment must be manufactured with milling equipment adequate for the designed geometry.
- Ensure proper fixation of the CORiTEC® preMill when it is being milled and ensure that the Interface of the CORiTEC® preMill is properly protected to prevent any damage to it. It is recommended that the CORiTEC® preMill Holder is used to achieve optimal fixing and proper protection.

Using digital workflow (intra oral scanning):

1. Import the digitized patient situation from the dental clinical intraoral scan to the design software by use of Elos Accurate® Scan Body in order to identify the position and orientation of the respective implant.
2. Import library file from Elos Medtech's webpage and select relevant implant platform from the library.
3. Design the customized abutment in the design software.
4. Export STL file to milling equipment.
5. Attach the CORiTEC® preMill to the milling equipment by use of the CORiTEC® preMill Holder and corresponding fixation screw.
6. Process and finalize the customized abutment in milling equipment according to the manufacturer's instruction for use.
7. Create a digital working model in the design software.
8. Export the STL file from the design software and send the STL file to your 3D printer or external 3D print provider.
9. Place an Elos Accurate® Model Analog for printed models in the 3D printed working model.

Using semi-digital workflow (desktop scanning):

1. Create working model and place an Elos Accurate® Model Analog appropriately in the model.
2. Place an Elos Accurate® Scan Body in the model analog.
3. By use of the Elos Accurate® Scan Body, scan the working model.
4. Remaining steps as for the digital workflow.

14.2 Creation and cementation of dental restoration (crown)

Recommended dental restoration material for CORiTEC® preMill (permanent restorations):

- Zirconia.

Compatible design software: 3Shape and Exocad.

Caution:

- Ensure to follow the selected dental restoration material minimum wall thickness and parameters according to the material supplier's instruction for use.
- The dental restoration must be manufactured with milling equipment adequate for the designed geometry.
- The dental restoration should if applicable be sintered in accordance with the instruction issued by the ceramic material supplier.

14.3 Preparing abutment and dental restoration for cementing

1. Blast bonding surface of the CORiTEC® preMill and restoration with aluminum oxide 50-150 µm and blasting pressure of 2 bar.
2. Clean the bonding surface of CORiTEC® preMill and restoration thoroughly with alcohol or with KATANA™ Cleaner by Kuraray Noritake.

14.4 Cementing the dental restoration

1. Seal the screw channels with wax.
2. Apply a primer to the bonding surface of CORiTEC® preMill and restoration and let it react according to manufactures instruction.

3. Apply self-adhesive dental cement on the CORiTEC® preMill . Follow the instruction for use of both the dental restoration material and cement material manufacturer. Bond the dental restoration to the CORiTEC® preMill placed in the working model or in the patients mouth. The milled dental restoration must fit the CORiTEC® preMill without gaps and voids.
4. Immediately remove excess cement from the os CORiTEC® preMill .

One of following bonding agents is recommended.

Primer / Cement	Restoration material	Intra oral cementation
Kuraray Noritake CLEARFIL™ CERAMIC PRIMER PLUS / PANAVIA™ V5	Zirconia	Yes
Ivoclar Vivadent® Monobond Plus / Multilink® Hybrid Abutment	Zirconia	No
3M™ Scotchbond™ Universal Plus Adhesive / RelyX™ Universal Resin Cement	Zirconia	Yes

Further instruction can be found in "Cementing guideline" at <https://elosdental.com/training-and-support/guides-our-products>

15 Use and handling for the dentist

The dentist receives the final dental restoration / working model with the Prosthetic screw from the dental lab.

1. Identify and unpack the Prosthetic screw(s).
2. Clean, disinfect and sterilize dental restoration and Prosthetic screw(s) according to this instruction for use.
3. Remove the healing cap, closure screw or temporary restoration from patient's mouth.
4. Gently, insert the dental restoration into the patient's mouth in proper position to the implant(s).
5. Place corresponding Prosthetic screw(s) in the dental restoration and tighten screws using torque values specified by the implant manufacturer.

16 Further information

For additional information about the use of Elos Medtech products, please contact your local sales representative.

17 Validity

Upon publication of this instruction for use, all previous versions are superseded.


18 Storage and Handling

CORiTEC® preMill should be stored at room temperature.

19 MRI Safety Information

 MRI Safety Information A person with an CORiTEC® preMill.Abutment, associated dental implant and prosthetic screw may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
Device Name	CORiTEC® preMill.
Static Magnetic Field Strength (B ₀)	1.5 T or 3.0 T
Maximum Spatial Field Gradient	20 T/m (2,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil
Operation Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR Image Artifact	The presence of passive implant devices may produce an image artifact that scales with the device size

	Do not use if package is damaged
	Non sterile
	Recommended torque
	Date of manufacture
	Medical Device
	Unique Device Identifier
	Keep dry
	Keep away from sunlight
	MR Conditional

20 Disposal

The dental restoration must be disposed as biological waste.

21 Serious Incidents:

If any serious incidents occur in relation with this device they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the country where the incident occurred.

A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat.

Contact Elos Medtech Pinol A/S at:

complaint.empi@elosmedtech.com

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. URL: <https://ec.europa.eu/tools/eudamed>.

22 Symbols

	Catalogue number
	Batch code
	Manufacturer
	Consult instructions for use
	Do not re-use
	Prescription only

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Şu adres üzerinden Elos Medtech Pinol A/S ile iletişime geçin: